



COMPANY PROFILES

제1공장 12726 경기도 광주시 곤지암읍 만삼로 108번길 63 (유사리 107) Tel. (031)763-8085, 8098 / 798-1136 Fax. (031)764-3526
제2공장 12726 경기도 광주시 곤지암읍 만삼로 80번길 (유사리 439-1) Tel. (031)763-8085, 8098 / 798-1136 Fax. (031)764-3526
제3공장 12726 경기도 광주시 곤지암읍 만삼로 108번길 (유사리 107) Tel. (031)763-8085, 8098 / 798-1136 Fax. (031)764-3526

http://www.kosanmedico.co.kr | http://www.kosanmedico.net | http://blog.naver.com/kosanco | E-mail: kosanco@hanmail.net



Best Business Partner

의료장비 전문제조업체

KOSANMEDICO



01 (주)고산메디코 (제1공장)



02 (주)고산메디코 (제2공장)



03 (주)고산메디코 (제3공장)



<http://www.kosanmedico.co.kr>

안녕하십니까 ?

눈부시게 발전하는 과학의 시대에서,

우리는 의료와 복지의 발달로 인하여, 건강하고 풍요롭게 살아가고 있습니다.

의료는 건강한 우리 생활에 반드시 필요한 부분이며, 건강이 높아지고 연령이 길어짐에 따라,

의료는 우리의 삶을 결정하는 매우 중요한 위치에 서게 되었습니다.

수요가 높아짐에 따라 준비된 의료 서비스 환경이 많이 필요로 하게 됩니다.

의료장비와 여러가지의 비품들은 효율적 의료 환경에 매우 필요한 요소들 입니다.

효율적인 환경을 구축하려면 잘 만들어야 하고, 잘 만들기 위해서는 경험과 기술이 뒷받침 되어야 합니다.

본 고산메디코는 수십년에 걸친 기술의 축적과 경험 등을 바탕으로,

여러분들이 필요로 하는 의료환경의 구축에 도움을 드리려 항상 노력합니다.

고객의 입장에서 필요한 전용 디자인 카트, 다품종 소량이 가능한 납품으로,

중소 사업자님들께 폭넓게 도움을 드리고저 합니다.

빠른 제작과 납기를 통해 수월한 사업력에 도움을 드리고, 추후 신속한 A/S까지 신용과 최선을 쌓도록 하겠습니다.

현재는 해외 시장을 선점하고자 한 걸음씩 나아가고 있습니다.

ISO 인증 등을 통하여 국내외에서 신뢰와 믿음을 보여드리려는 그런 기업이 되도록 노력하고 있습니다.

앞으로, 해외 어느 곳에서든지 본 고산메디코의 제품들을 보실 수 있을 것 입니다.

앞으로도 더 많은 노력을 기울여서 의료산업 분야의 선두주자가 되어 여러분의 성공 파트너로써,

좋은품질, 믿을 수 있는 기업으로 더 높게 성장하도록 하겠습니다.

여러분들의 든든한 파트너가 되어 드리겠습니다.

감사합니다.

-주식회사 고산메디코 대표 이상철-


연혁

| | |
|-----------|-----------------------------|
| 1994년 8월 | 개인사업자 설립 |
| 2000년 2월 | 현 소재지 신축이전 |
| 2006년 1월 | 법인전환 |
| 2009년 5월 | 벤처기업 확인 <기술보증기금> |
| 2011년 3월 | iso9001:2008 <ICR> |
| 2011년 9월 | 산학교류협약 <한국폴리텍 I대학 성남캠퍼스> |
| 2012년 2월 | KIMES 출품 |
| 2012년 6월 | 의료용카트 특허 <특허청> |
| 2012년 9월 | 한국무역협회 회원 <한국무역협회> |
| 2012년 10월 | ENDOSCOPE350 CE인증 <TUVNORD> |
| 2013년 3월 | KIMES 출품 |
| 2014년 3월 | KIMES 출품 |
| 2015년 3월 | KIMES 출품 |
| 2015년 10월 | 2공장 설립 |
| 2017년 12월 | 3공장 설립 |

서울의료원, 서울성모병원, 충남대학교병원, 국제성모병원 등 국내 중대형 병원 납품



- 기존의 모델을 구매 및 사용하실 수 있습니다.
- 기존의 모델을 수정 할 수 있습니다.
- 전용 모델을 제작하실 수 있습니다.



KOSANMEDICO

| | |
|---------|---|
| 회사명 | 주식회사 고산메디코 (KOSANMEDICO.CO.,LTD) |
| 사업자등록번호 | 126-81-91869 |
| 연 혁 | 1994년 8월 1일 설립 2000년 2월 1일 신축이전 2006년 1월 2일 법인전환 |
| 대표자 | 이상철 |
| 전 화 | 031-763-8085, 031-763-8098, 031-798-1136 |
| 팩 스 | 031-764-3526 |
| e-mail | kosanco@hanmail.net |
| 웹사이트 | http://www.kosanmedico.co.kr , http://www.kosanmedico.net |
| 주 소 | 464-871 경기도 광주시 곤지암읍 만삼로 108번길 63 |

KOSANMEDICO





Contents Page

수술실비품

Operating Room Appliances

| | | | |
|---------------------------------|-------|--------------------|----|
| INSTRUMENT CABINET | 8-13 | WASH BASIN STAND | 21 |
| SCRUB STATION | 14-15 | RACK DRYING GLOVE | 21 |
| SINK CABINET | 15 | SUS WIRE BASKET | 21 |
| WALL-MOUNTED INSTRUMENT CABINET | 16-18 | SUCTION BOTTLE | 22 |
| PACS CABINET | 18 | CHEST BOTTLE STAND | 22 |
| RACK BASKET SYSTEM | 18 | KICK BUCKET | 23 |
| RACK SYSTEM | 18 | TUBE STAND | 24 |
| INSTRUMENT TABLE | 19 | SIDE LAMP | 24 |
| MAYO STAND | 20 | MONITOR CART(SUS) | 24 |
| MAYO TRAY | 20 | | |

간호실비품

Nursing Appliances

| | | | |
|-----------------------------|-------|------------------------------------|-------|
| MEDICAL CART | 26-27 | DRESSING CART | 51 |
| EMERGENCY CART | 28-31 | DRESSING CART (BANDING TYPE) | 51 |
| MEDICATION CART (DRAWER) | 32-37 | TREATMENT CART | 52 |
| MEDICATION CART (CASSETTE) | 38-44 | TREATMENT CART (BANDING TYPE) | 52 |
| CASSETTE DELIVERY CART | 44 | SIDE TABLE | 53 |
| BIN BOX | 45 | BED PEN RACK | 53 |
| ANESTHESIA CART | 46-47 | URINE CONTAINER RACK (HOLDER TYPE) | 53 |
| ICU CART | 47 | FECES PEN RACK (HOLDER TYPE) | 53 |
| SUPER PLUS CART | 47 | LINEN HAMPER CART | 54 |
| UTILITY CART | 47 | LINEN HAMPER BAG | 54 |
| HAIR SHAMPOOING MACHINE | 48 | STEP FOOT STOOL | 55 |
| WALKER | 48 | CHART HOLDER CART | 55 |
| WHEEL CHAIR | 48 | BASKET CART | 56-57 |
| HAIR SHAMPOOING BASKET | 48 | OXYGEN TANK HOLDER | 56 |
| UTILITY CART (DRAWER TYPE) | 49 | OXYGEN TANK CART | 57 |
| UTILITY CART (A,B,S DRAWER) | 50 | SCREEN | 58 |
| UTILITY CART | 51 | 환자운반기(유압식) | 59 |
| | | 환자운반기(TOTAL LOCK) | 59 |

진료실비품

Consultation Room Appliances

| | | | |
|---------------------------------|----|-----------------------|-------|
| ENDOSCOPIC CABINET | 60 | ENDOSCOPIC LENS SAFER | 64 |
| COLONOFIBERSCOPE CABINET | 60 | APRON STORAGE RACK | 64 |
| ENDOSCOPIC CABINET(ROTATE TYPE) | 61 | X-RAY ILLUMINATOR | 64 |
| ACCESSORIES STAND | 62 | FILM CART | 65 |
| ENDOSCOPIC DRYING RACK STAND | 62 | OPERATING STOOL | 66-67 |
| ENDOSCOPIC DRYING RACK (WOOD) | 63 | | |

공급실비품

Supply Room Appliances

| | | | |
|---------------------------|-------|-------------------------|-------|
| WIRE SHELF CART | 68 | DUST BOX CART | 74 |
| WIRE SHELF | 69 | DELIVERY CART | 74 |
| PLASTIC DOUBLE CART | 70 | FOOD CART | 74 |
| PLASTIC TRIPLE CART | 70 | SAFETY MEDICINE CABINET | 75 |
| SUDDING METAL RACK SYSTEM | 71 | MULTI BIN(PLASTIC) | 76-80 |
| TRUCK CART | 72 | MULTI BIN CART | 81 |
| DANPRA(OPTION) | 72 | MULTI BIN CABINET(WOOD) | 82-83 |
| LINEN SERVICE CART | 72,74 | MULTI BIN CART | 84 |
| WORK TABLE | 73 | TRAY RACK CART (SUS) | 84 |
| SEALING TABLE | 74 | PASS BOX SYSTEM | 84 |
| SEALING STAND | 74 | U,V STERILIZER | 84 |

의료용침대

Hospital Bed

| | | | |
|-------------------------------|-------|---------------------------|----|
| BED MONITOR ARM | 85 | EXAMINING BED | 94 |
| MOTORIZED PATIENT-CONTROL BED | 86-89 | SHOWER TROLLEY (STANDARD) | 94 |
| HOSPITAL BED | 90 | INFANT BASSINET STAND | 95 |
| 쇼아용 침대 | 91-92 | STRETCHER CART | 95 |
| PILLOW/MATRESS | 93 | 3단 STRETCHER CART | 96 |
| ATTENDANT BED | 93 | 3단 PROTECTOR BED/CHAIR | 96 |
| EXAMINING TABLE | 93 | SHOWER TROLLEY | 97 |
| CABINET (WOOD) | 94 | | |

카트

Carts

| | | | |
|---------------------|---------|-----------------------------|---------|
| ENDOSCOPE CART | 98-101 | 혈압계용 | 111 |
| ASSIST MONITOR CART | 98 | EKG CART | 112 |
| ARTHROSCOPE CART | 101-102 | STANDARD HANGER | 112-114 |
| LAPAROSCOPE CART | 102-104 | JOINT HANGER | 112 |
| E,N,T CART | 104-105 | NEBULIZER | 114 |
| PC CART | 105-108 | BEAUTY | 114 |
| MULTI CART | 108-110 | MEB | 114 |
| BOVIE CART | 111 | ASSIST MONITOR CART (MOTOR) | 115 |

모니터카트

Monitor Carts

| | |
|---------------|---------|
| Monitor Carts | 116-119 |
|---------------|---------|

수액걸이대

IV Pole Stand

| | |
|---------------|---------|
| IV Pole Stand | 120-123 |
|---------------|---------|

모니터암

Monitor Arm

| | | | |
|---------------|-----|----------------------|-----|
| S-TYPE | 124 | EXTENSION HOLDER | 128 |
| I-TYPE | 125 | CENTER HEAD (1 HOLE) | 128 |
| JOINT ARM | 126 | CENTER HEAD (2 HOLE) | 128 |
| EXTENSION ARM | 127 | DUAL MONITOR ARM | 128 |
| CENTER ARM | 127 | STAND MONITOR CARM | 128 |
| EXTENSION 160 | 128 | TABLET HOLDER | 129 |
| EXTENSION 250 | 128 | NOTEBOOK PLATE | 129 |
| EXTENSION 450 | 128 | SINGLE MONITOR ARM | 129 |

월마운트

Wall Mount

| | | | |
|---------------------|---------|-------------|-----|
| WALL MOUNT | 130-131 | TABLE MOUNT | 131 |
| WALL MOUNT (BASKET) | 130 | | |

폴마운트

Pole Mount

| | | | |
|------------|-----|-------------|-----|
| POLE MOUNT | 132 | PLATE MOUNT | 132 |
|------------|-----|-------------|-----|

악세서리 - 바구니,열티랩,기타

Accessory

| | | | |
|-------------|-----|-----------|---------|
| WIRE BASKET | 133 | MULTI-TAP | 133-134 |
| MASH BASKET | 133 | | |

A.B.S Cart

A.B.S Cart

| | |
|------------|---------|
| A.B.S Cart | 138-161 |
|------------|---------|

주문형 카트

Order Type Carts

| | |
|------------------|---------|
| Order Type Carts | 162-173 |
|------------------|---------|

다양한 규격의 제품도 제작이 가능함 (주문제작)

- Free-standing instrument cabinets are designed for larger storage space and Constructed of heavy-duty stainless steel.
- Body Frame : All Stainless Steel 304 1.0t & 1.2t Plate Production.
- Glass : 5.0t Strong Glass.
- Shelf : Two shelves is adjustable in 30mm intervals.
- Handle : Locking Key unit.
- Adjustable Floor Level is Possible



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|-----------|
| KS-IC3 | 900×450×1800 | 여닫이(Door) |

* Four Shelves is Adjustable in 30mm Intervals



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|-----------|
| KS-IC4 | 900×450×1800 | Open type |
| KS-IC4-1 | 1200×450×1800 | Open type |



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-IC1 | 1200×500×1800 | 주문제작가능 |

* Bin : KS-SMB504 12EA



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-IC2 | 1200×500×1800 | 주문제작가능 |



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|-----------|
| KS-IC5 | 900×450×1800 | 여닫이(Door) |

* Four Shelves is Adjustable in 30mm Intervals



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|-----------|
| KS-IC6 | 1200×500×1800 | 미닫이(Door) |

* Four Shelves is Adjustable in 30mm Intervals



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|-----------|
| KS-IC7 | 1200×500×1800 | 미닫이(Door) |

* Four Shelves is Adjustable in 30mm Intervals



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|-----------|
| KS-IC8 | 900×450×1800 | 여닫이(Door) |

* Four Shelves is Adjustable in 30mm Intervals



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|-----------|
| KS-IC11 | 900×450×1800 | 미닫이(Door) |
| KS-IC11-1 | 1200×500×1800 | 미닫이(Door) |



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|-----------|
| KS-IC12 | 900×450×1800 | 여닫이(Door) |



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|-----------|
| KS-IC9 | 900×450×1800 | 여닫이(Door) |

* Two Shelves is Adjustable in 30mm Intervals
* Multi Plastic Bin : 18EA



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|-----------|
| KS-IC10 | 900×450×1800 | 여닫이(Door) |

* Three Shelves is Adjustable in 30mm Intervals
* Drawer : Slide Ball Bearing Rail



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-IC13 | 1200×500×1800 | - |



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-IC14 | 1200×500×1800 | - |



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-IC15 | 1300×550×1800 | — |


■ Composition Units : SUS Shelf (2EA)
ABS Basket (8EA), SUS Wire Basket (2EA)



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-IC16 | 1300×550×1800 | — |

■ Composition Units : SUS Shelf (2EA)
ABS Basket (10EA), SUS Wire Basket (4EA)



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-IC19 | 900×450×1800 | Open |

■ Composition Units : SUS Shelf (3EA)
Multi Plastic Bin (30EA)



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-IC20 | 1300×550×1800 | — |


■ Composition Units : Multi Plastic Bin (24EA)
ABS Basket (6EA), SUS Wire Basket (4EA)



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|-----------|
| KS-IC17 | 1300×550×1800 | 미닫이(Door) |


■ Composition Units : ABS Basket (16EA), SUS Wire Basket (8EA)



INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-IC18 | 700×550×1800 | OPEN |

■ Composition Units : ABS Basket (8EA), SUS Wire Basket (4EA)




INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-IC21 | 1200×500×1800 | 보조판 있음 |



INSTRUMENT CABINET


| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-IC22 | 1200×500×1800 | 보조판 없음 |



SCRUB STATION

| Model | Size (L×D×H / mm) | Soap Pump |
|------------|-------------------|------------|
| KS-PMS-100 | 700×600×1150 | 1(220V-3W) |
| KS-PMS-200 | 1400×600×1150 | 1(220V-3W) |
| KS-PMS-300 | 2000×600×1150 | 2(220V-3W) |

- Treatment : Germs removal by filtration
- Heater : 2.0Kw×1EA
- Water : 20~60℃
- Water time : Adjustable Within 30sec
- Filter : Pre-Filter 1EA, Main-Filter 1EA
- Sterilizer : Ultra Violet Ray Flow-Sterilizing System (15W x 2EA)
- U.V Lamp : 2,000 Hours, Using Time indicates on control panel with digital (Filter life time check)
- Water supply : Min 1.0kg/㎥ Max 5.0kg/㎥
- Water temperature control : Thermostatic Mixing Valve
- Soap filling check : Microprocessor Auto Control
- Power source : AC 220V, 50/60Hz
- Power consumption : 2.2Kw
- Water tank capacity : Approx 30L
- Net weight [kg] 75 / 105 / 135



KS-PMS-100

KS-PMS-200

KS-PMS-300

SPECIAL FEATURE



1. Attach UV Lamp check LED Lamp to the front of the equipment



2. Splash bump design



3. Install pollution control soap case

APPLICATION

- Surgical operating room
- Delivery room
- Emergency room
- ICU room
- Radiology room
- Other places in need of sterilized water

FEATURE

- Heavy duty stainless steel 304 construction for rust proof and corrosion resistance
- Sterilization system by UV lamps and dual filterings by main filter 5um and pre-filter 10um
- Sterilized water of over 99% is provided
- Easy operation by knee-water and soap
- Water pressure control device

BASIC OPERATION AND FUNCTION

1.Heating

- 1) Hot water supply in winter
- 2) Fast warming by one heater (2KW)
- 3) Two safety thermo-switches to control the in-creasin pressure of the warming tank
- 4) Constant water temperature by safety thermo-switches

2.Sterilization

- 1) The water can be sterilized to 99.5%(0.6T/H) by UV water sterilizing lamps
- 2) The UV lamps are made of dual quart tubes for electrical safety
- 3) Lamp us time can be read as life-timer is attached

3.Filtering

- 1) pre-filter 10um
- 2) Main-filter 5um
- 3) Thorough filtering by dual-filters
- 4)Easy filter replacement



SINK CABINET

| Model | Size (L×D×H / mm) | Remark |
|---------|--------------------|--------|
| KS-SK01 | 1200×600×(850+100) | - |
| KS-SK02 | 1500×600×(850+100) | - |
| KS-SK03 | 1800×600×(850+100) | - |

■ SINK CABINET 주문 제작 가능합니다.

SINK

넓은 면적의 내식성이 우수한 스테인리스 싱크는 각종 의료기구를 편리하게 세척 할 수 있습니다.



WALL-MOUNTED INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|-------------|
| KS-WMC01 | 3700×500×2300 | Top Cabinet |

■ 상기 타입 이외에도 병원 실정에 적합하게 주문 제작 및 설치가 가능

COMPOSITION UNITS

- ① Basket Cabinet Accessories
 - ABS Basket (Low Type) 2EA
 - ABS Basket (High Type) 4EA
 - SUS Wirt Basket (Low Type) 2EA
 - SUS Wirt Basket (High Type) 2EA
 - Multi Box 24EA (Partition)
 - Supporting Slide Table
- ② Instrument Cabinet Accessories,
 - Shelf 2EA
- ③ Instrument Cabinet Accessories,
 - Shelf 2EA
 - Multi Plastic Bin 24EA (Partition)
- ④ Supporting Slide Table
- ⑤ The Top Cabinet



WALL-MOUNTED INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|-------------|
| KS-WMC03 | 3400×500×2300 | Top Cabinet |



WALL-MOUNTED INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|-------------|
| KS-WMC04 | 3200×500×2300 | Top Cabinet |



WALL-MOUNTED INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|-------------|
| KS-WMC02 | 4500×500×2300 | Top Cabinet |

■ 상기 타입 이외에도 병원 실정에 적합하게 주문 제작 및 설치가 가능

COMPOSITION UNITS

- ①② Instrument Drawer Cabinet
 - Shelf 2EA
- ③ Instrument Bin Cabinet Accessories,
 - Multi Box (Low) 34EA (Partition)
 - Multi Box (High) 32EA (Partition)
- ④ Instrument Cabinet Accessories,
 - Shelf 2EA
- ⑤ The Top Cabinet
- ⑦ Drawer (Swivels Sliding Rail)



WALL-MOUNTED INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-WMC05 | 4200×500×1800 | - |



WALL-MOUNTED INSTRUMENT CABINET

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|-------------|
| KS-WMC06 | 2400×500×2300 | Top Cabinet |

■ 상기 타입 이외에도 병원 실정에 적합하게 주문 제작 및 설치가 가능

COMPOSITION UNITS

- ① Basket Cabinet Accessories
 - ABS Basket (Low Type) 2EA
 - ABS Basket (High Type) 4EA
 - SUS Wire Basket (Low Type) 2EA
 - SUS Wire Basket (High Type) 2EA
 - Multi Box 24EA (Partition)
 - Supporting Slide Table.
 ② Instrument Cabinet Accessories.
 - Shelf 2EA
 ③ The Top Cabinet



RACK BASKET SYSTEM

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|-----------------------------------|
| KS-RBS | 2450×550×1800 | STS Wire Basket 560×410×125(H) |



PACS CABINET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|-------------|
| KS-PC01 | 1200×500×2300 | Top Cabinet |

■ Option : Monitor, Viewbox (2Panel)



RACK SYSTEM

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-RACK100 | 추후 협의 | |



INSTRUMENT TABLE (ROUND TYPE)

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-ITR1 | 1200×500×900 | Round |
| KS-ITR2 | 1600×500×900 | Round |
| KS-ITR3 | 1900×500×900 | Round |

■ All Stainless Steel 304 Construction
 ■ 4" Swivel Brake Caster 6EA
 ■ Two Self shaped with guard rail



INSTRUMENT TABLE (V-TYPE)

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-ITV1 | 1200×500×900 | V형 |
| KS-ITV2 | 1600×500×900 | V형 |
| KS-ITV3 | 1900×500×900 | V형 |

■ All Stainless Steel 304 Construction
 ■ 4" Swivel Brake Caster 6EA
 ■ Two Self shaped with guard rail



INSTRUMENT TABLE (I-TYPE)

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-ITI1 | 900×500×900 | 일자형 |
| KS-ITI2 | 1200×500×900 | 일자형 |
| KS-ITI3 | 1600×500×900 | 일자형 |

■ All Stainless Steel 304 Construction
 ■ 4" Swivel Brake Caster 4EA (Brake 2EA)
 ■ 상기 모델은 주문에 따라 규격변동이 가능하며 상부 선반에 GUARD RAIL 제작 가능함



* O.R (수술실용)

* 외래 진료실용

KS-MSA-1 600×400×8
= O.R (수술실용) Foot Switch(자동)



- All Stainless Steel Pipe Frame
- Stainless Basin
- 2" Castor 4EA



- All Stainless Steel Pipe Frame
- Stainless Basin
- 2" Caster 4EA



* 외래 진료실용



| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-RDG20 | 600×400×1500 | 20高 |
| KS-RDG30 | 800×400×1500 | 30高 |
| KS-RDG50 | 1200×400×1500 | 50高 |



| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|---------|
| KS-SWB-1 | 560×410×130 | SUS Bar |
| KS-SWB-2 | 560×410×190 | SUS Bar |



SUCTION BOTTLE 3000 CC

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-SB-30C | φ350×500 | Cap 포함 |

- * Plastic 3000cc bottle (Polycarbonate)
- * Sus bottle cap : 2 air line
- * Option accessories : Silicon horse
- * Overflow function

SUCTION BOTTLE 8000 CC

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-SB-80C | φ350×550 | Bottle |
| KS-SB-80CA | — | Cap |

- * Plastic 8000cc bottle (Polycarbonate)
- * Sus bottle cap : 2 air line
- * Option accessories : Silicon horse
- * Overflow function

SUCTION BOTTLE HOLDER

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-SB-H80C | 600×350×500 | — |

- * Swivel caster 2" 4EA
- * Urethane bumper
- * Plastic 8000cc bottle (Polycarbonate)
- * Sus bottle cap : 2 air line
- * Silicon horse (2M-2EA)



CHEST BOTTLE STAND

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-LBS-01 | 180×100×285 | — |

- * All sus304 construction



CHEST BOTTLE STAND

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-LBS-02 | 165×145×285 | — |

- * All sus304 construction



KICK BUCKET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-KB02 | φ350×350 | — |

- * Accessories Stainless Steel Bucket



KICK BUCKET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|-------------|
| KS-KB03 | φ350×700 | Handle Type |


- * Accessories Stainless Steel Bucket



KICK BUCKET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-KB05 | φ430×900 | — |

- * Accessories Stainless Steel Bucket



KICK BUCKET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-KB04 | φ350×140 | — |

- * Accessories Steel Bucket



TUBE STAND

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-TS1 | 210×170×500 | STS304 |

TUBE STAND

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-TS2 | 300×200×180 | STS/AL |

SIDE LAMP

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|------------|
| KS-SL02 | 1430×2000 | 25,000 LUX |



MONITOR CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-MO1 | 600×600×1300 | - |

MONITOR CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-MO2 | 550×500×1300 | - |

MONITOR CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-MO3 | 600×500×1200 | - |

■ MONITOR CART는 주문에 따른 제작이 가능함





환자 상태의 기록을 위한 다목적용 카트입니다. 노트북 및 데스크탑 PC의 설치 및 사용이 가능합니다. 데스크탑 PC의 사용시, 별도의 모니터 암(ARM)의 사용이 가능합니다. 수납칸의 자유로운 칸 배열 가능합니다. 환경에 따라 필요한 구성품들을 추가할 수 있습니다.



MEDICAL CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-MDC-D1 | 740×575×900/1000 | - |
| KS-MDCS-D1 | 610×575×900/1000 | - |

* OPTION 높이

MEDICAL CART

| Model | Size (L×D×H / mm) | Remark |
|---------------|-------------------|--------|
| KS-MDC-D1-P1 | 740×575×900/1000 | - |
| KS-MDC-D1S-P1 | 610×575×900/1000 | - |

적용 가능 구성품 (Option)



MONITOR ARM

| Model | Remark |
|------------|----------|
| KS-ARM-MSB | A,B,S 전용 |

* HEAD 포함
* OPTION



BASKET

| Model | Size (L×D×H / mm) |
|---------|-------------------|
| KS-AB-1 | 560×410×55 |
| KS-AB-2 | 560×410×125 |
| KS-AB-3 | 560×410×190 |

* DIVIDER를 통한 칸 조절
* OPTION



IV POLE

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-IV05 | Φ22 |

* H=450
* 5 HOOKS
* OPTION



NEEDLE HOLDER

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-NH02 | 155×120(110) |


* OPTION



WASTE BOX

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-WB01 | 185×175×330 |

* OPTION



WASTE HOLDER

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-WH02 | 210×230×155 |

* 신형모델 적용
* OPTION



MEDICAL CART

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-MDC | 740×575×900/1000 | - |
| KS-MDCS | 610×575×900/1000 | - |



MEDICAL CART

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-MDC-YD1 | 740×575×900/1000 | - |
| KS-MDCS-YD1 | 610×575×900/1000 | - |



MEDICAL CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-MDC-P1 | 740×575×900/1000 | - |
| KS-MDCS-P1 | 610×575×900/1000 | - |



MEDICAL CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-MDC-D2 | 740×575×900/1000 | - |
| KS-MDCS-D2 | 610×575×900/1000 | - |



MEDICAL CART

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-MDCS-D2W2 | 740×575×900/1000 | - |



MEDICAL CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-MDC-D3 | 740×575×1000 | - |
| KS-MDCS-D3 | 610×575×1000 | - |



MEDICAL CART

| Model | Size (L×D×H / mm) | Remark |
|----------------|-------------------|----------|
| KS-MDC-KD1-P1 | 740×575×900/1000 | Keyboard |
| KS-MDCS-KD1-P1 | 610×575×900/1000 | Keyboard |



MEDICAL CART

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|----------|
| KS-MDC-KD2 | 740×575×900/1000 | Keyboard |
| KS-MDCS-KD2 | 610×575×900/1000 | Keyboard |



IV POLE


| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-AC-IV06 | Φ22 | - |



응급환자의 상태 관리를 위한 다목적용 카트입니다. 제세동기 장치의 거치가 가능하며 약의 종류별 보관에 용이합니다. 메모 기록을 위한 별도의 슬라이드 판이 있으며 다양한 구성품을 적용하여 보다 환경에 맞게 사용할 수 있습니다.

EMERGENCY CART

- Size : 740×575×1250(2200)(L×D×H / mm)
- LV Pole 5 hooks
- Defibrillator Plate
- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- C.P.R Board



서랍을 열고 바스켓을 꺼낼 수 있습니다. 보관된 내용물과 함께 이동이 가능합니다.


OXYGEN TANK

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-OT01 | For Φ120 |

Oxygen tank 와 C.P.R Board의 장착이 가능합니다. (옵션사항)

C.P.R BOARD

| Model | Size (L×D×H / mm) |
|-------------|-------------------|
| KS-AC-CPR01 | 420×530×20 |



공통적용 가능 구성품 (OPTION)



DEFIBRILLATOR


| Model | Size (L×D×H / mm) |
|------------|----------------------------|
| KS-AC-DF01 | L : 395×275 S : 330×235 |

= OPTION



SLIDE PLATE

| Model | Size (L×D×H / mm) |
|-------|-------------------|
| - | 350×340 |



BASKET

| Model | Size (L×D×H / mm) |
|---------|-------------------|
| KS-AB-1 | 560×410×55 |
| KS-AB-2 | 560×410×125 |
| KS-AB-3 | 560×410×190 |

= DIVIDER를 통한 칸 조절
= OPTION



UTILITY CONTAINER

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-UC01 | 340×100×300(250) |

= OPTION



STORAGE BOX

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-SB01 | 480×180×418 |

= 6 PERSON
= OPTION
= KS-CA703



IV POLE

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-IV05 | Φ22 |

= H=450
= 5 HOOKS
= OPTION



NEEDLE HOLDER

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-NH02 | 165×130×110 |

= OPTION
= 2L



WASTE BOX

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-WB01 | 185×175×330 |

= OPTION



WASTE HOLDER

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-WH02 | 210×230×155(5L) |
| KS-AC-WH03 | 280×300×155(12L) |

= 신형모델 적용
= OPTION



EMERGENCY CART (740)

| Model | Size (L×D×H / mm) | Remark |
|----------|--------------------|----------|
| KS-EMC5D | 740×575×1000(1070) | Standard |

※ Lock (Key)



EMERGENCY CART (740)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|--------|
| KS-EMC5D-01 | 740×575×1000(1070) | — |

※ Utility Container
 ※ Waste Holder
 ※ Waste Box
 ※ Needle Holder
 ※ Lock (Key)



EMERGENCY CART (610)

| Model | Size (L×D×H / mm) | Remark |
|------------|--------------------|----------|
| KS-EMC5D-S | 610×575×1000(1070) | Standard |



EMERGENCY CART (610)

| Model | Size (L×D×H / mm) | Remark |
|--------------|--------------------|--------|
| KS-EMC5D-01S | 610×575×1000(1070) | — |

※ Utility Container
 ※ Waste Holder
 ※ Waste Box
 ※ Needle Holder



EMERGENCY CART (740)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|--------|
| KS-EMC5D-02 | 740×575×1000(1070) | — |

※ Storage Box (6 person)
 ※ Waste Holder
 ※ Waste Box
 ※ Needle Holder
 ※ Lock (Key)



EMERGENCY CART (740)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|--------|
| KS-EMC5D-03 | 740×575×1000(1070) | — |

※ Storage Box (6 person) 2ea
 ※ Lock (Key)



EMERGENCY CART (610)

| Model | Size (L×D×H / mm) | Remark |
|--------------|--------------------|--------|
| KS-EMC5D-02S | 610×575×1000(1070) | — |

※ Storage Box (6 person)
 ※ Waste Holder
 ※ Waste Box
 ※ Needle Holder



EMERGENCY CART (610)

| Model | Size (L×D×H / mm) | Remark |
|--------------|--------------------|--------|
| KS-EMC5D-03S | 610×575×1000(1070) | — |


※ Storage Box (6 person) 2ea



환자의 약물 치료를 위한 간호사용 카트입니다. 이동성과 수납을 위해 수납함을 앞뒤로 배치됩니다. 여러 종류의 칸을 이용하여 적절한 분배가 가능합니다. 한번의 잠금 조작으로 모든 칸을 잠금 가능합니다. 각 Bin Box는 수납의 조합이 가능합니다. 확장을 통하여 관리대상의 숫자만큼 확장 가능합니다.

MEDICATION CART (DRAWER)

- Size : 740×575×1250(2200)(L×D×H / mm)
- LV Pole 5 hooks
- Notebook Plate
- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Bin Box 30 Person



Bin Box의 내부는 칸막이를 이용하여 공간을 조절할 수 있습니다.



Key를 이용하여 Bin Box를 전체적으로 잠금하여 분실에 대비합니다.

공통적용 가능 구성품 (OPTION)



NOTEBOOK PLATE

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-NT01 | 395×275 |

= OPTION



SLIDE PLATE

| Model | Size (L×D×H / mm) |
|-------|-------------------|
| - | 350×340 |



UTILITY CONTAINER

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-UC01 | 340×100×300(250) |

= OPTION



STORAGE BOX

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-SB01 | 480×180×418 |

= 6 PERSON
= OPTION
= KS-CA703



IV POLE

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-IV05 | Φ22 |

= H=450
= 5 HOOKS
= OPTION



NEEDLE HOLDER

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-NH02 | 165×130×110 |

= OPTION
= 28



WASTE BOX

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-WB01 | 185×175×330 |

= OPTION



WASTE HOLDER

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-WH02 | 210×230×155 |

= 신형모델 적용
= OPTION



MEDICATION CART (740,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|----------|--------------------|-----------|
| KS-MCC20 | 740×575×1000(1070) | 20 Person |

- Standard
- Lock (Key)

MEDICATION CART (740,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|----------|--------------------|-----------|
| KS-MCC30 | 740×575×1000(1070) | 30 Person |

- Standard
- Lock (Key)

MEDICATION CART (740,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|----------|--------------------|-----------|
| KS-MCC40 | 740×575×1000(1070) | 40 Person |

- Standard
- Lock (Key)



MEDICATION CART (740,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC20-02 | 740×575×1000(1070) | 20 Person |

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)

MEDICATION CART (740,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC30-02 | 740×575×1000(1070) | 30 Person |

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)

MEDICATION CART (740,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC40-02 | 740×575×1000(1070) | 40 Person |

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)



MEDICATION CART (740,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC20-03 | 740×575×1000(1070) | 20 Person |

- Storage Box (6 person) 2ea
- Lock (Key)

MEDICATION CART (740,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC30-03 | 740×575×1000(1070) | 30 Person |

- Storage Box (6 person) 2ea
- Lock (Key)

MEDICATION CART (740,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC40-03 | 740×575×1000(1070) | 40 Person |

- Storage Box (6 person) 2ea
- Lock (Key)



MEDICATION CART (740,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC20-01 | 740×575×1000(1070) | 20 Person |

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)

MEDICATION CART (740,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC30-01 | 740×575×1000(1070) | 30 Person |

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)

MEDICATION CART (740,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC40-01 | 740×575×1000(1070) | 40 Person |

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)



MEDICATION CART (610,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|----------|--------------------|-----------|
| KS-MCC16 | 610×575×1000(1070) | 16 Person |

- Standard
- Lock (Key)



MEDICATION CART (610,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|----------|--------------------|-----------|
| KS-MCC24 | 610×575×1000(1070) | 24 Person |

- Standard
- Lock (Key)



MEDICATION CART (610,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|----------|--------------------|-----------|
| KS-MCC32 | 610×575×1000(1070) | 32 Person |

- Standard
- Lock (Key)



MEDICATION CART (610,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC16-02 | 610×575×1000(1070) | 16 Person |

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)



MEDICATION CART (610,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC24-02 | 610×575×1000(1070) | 24 Person |

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)



MEDICATION CART (610,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC32-02 | 610×575×1000(1070) | 32 Person |

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)



MEDICATION CART (610,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC16-01 | 610×575×1000(1070) | 12 Person |

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)



MEDICATION CART (610,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC24-01 | 610×575×1000(1070) | 24 Person |

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)



MEDICATION CART (610,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC32-01 | 610×575×1000(1070) | 32 Person |

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)



MEDICATION CART (610,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC16-03 | 610×575×1000(1070) | 16 Person |

- Storage Box (6 person) 2ea
- Lock (Key)



MEDICATION CART (610,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC24-03 | 610×575×1000(1070) | 24 Person |

- Storage Box (6 person) 2ea
- Lock (Key)



MEDICATION CART (610,DRAWER)

| Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|-----------|
| KS-MCC32-03 | 610×575×1000(1070) | 32 Person |



- Storage Box (6 person) 2ea
- Lock (Key)



환자의 약물 치료를 위한 간호사용 카트입니다. 이동성을 위해 카세트 방식으로 제작되어 카트의 이동 이외에 단독으로 Bin Box의 이동 및 활용이 가능합니다. 한번의 잠금 조작으로 모든 칸을 잠금 가능합니다. 후면에도 Multi Bin이 존재하여 수납의 기능을 하고 있습니다.

MEDICATION CART (CASSETTE)

- Size : 740×575×1250(2200) (L×D×H / mm)
- LV Pole 5 hooks
- Notebook Plate
- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Bin Box 10 Person
- Cassette Box
- Multi Bin 10 Person



공통적용 가능 구성품 (OPTION)



NOTEBOOK PLATE


| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-NT01 | 395X275 |

= OPTION



SLIDE PLATE

| Model | Size (L×D×H / mm) |
|-------|-------------------|
| - | 350X340 |



BASKET

| Model | Size (L×D×H / mm) |
|---------|-------------------|
| KS-AB-1 | 560×410×55 |
| KS-AB-2 | 560×410×125 |
| KS-AB-3 | 560×410×190 |

= DIVIDER를 통한 칸 조절
= OPTION



UTILITY CONTAINER

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-UC01 | 340X100X300(250) |

= OPTION



STORAGE BOX

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-SB01 | 480X180X418 |

= 6 PERSON
= OPTION
= KS-CA703



IV POLE

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-IV05 | Φ22 |

= H=450
= 5 HOOKS
= OPTION



NEEDLE HOLDER

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-NH02 | 165×130×110 |

= OPTION
= 20



WASTE BOX

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-WB01 | 185×175×330 |

= OPTION



WASTE HOLDER

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-AC-WH02 | 210×230×155 |

= 신형모델 적용
= OPTION



MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE)

| Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark |
|-----------|--------------------|-----------|-----------|--------------------|-----------|-----------|--------------------|-----------|
| KS-MCC10C | 740×575×1000(1070) | 10 Person | KS-MCC15C | 740×575×1000(1070) | 15 Person | KS-MCC20C | 740×575×1000(1070) | 20 Person |

- Standard
- Lock (Key)
- Cassette Type

- Standard
- Lock (Key)
- Cassette Type

- Standard
- Lock (Key)
- Cassette Type



MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE)

| Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark |
|--------------|--------------------|-----------|--------------|--------------------|-----------|--------------|--------------------|-----------|
| KS-MCC10C-01 | 740×575×1000(1070) | 10 Person | KS-MCC15C-01 | 740×575×1000(1070) | 15 Person | KS-MCC20C-01 | 740×575×1000(1070) | 20 Person |

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)
- Cassette Type

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)
- Cassette Type

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)
- Cassette Type



MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE)

| Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark |
|--------------|--------------------|-----------|--------------|--------------------|-----------|--------------|--------------------|-----------|
| KS-MCC10C-02 | 740×575×1000(1070) | 10 Person | KS-MCC15C-02 | 740×575×1000(1070) | 15 Person | KS-MCC20C-02 | 740×575×1000(1070) | 20 Person |

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)
- Cassette Type

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)
- Cassette Type

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)
- Cassette Type



MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE) MEDICATION CART (740,CASSETTE)

| Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark |
|--------------|--------------------|-----------|--------------|--------------------|-----------|--------------|--------------------|-----------|
| KS-MCC10C-03 | 740×575×1000(1070) | 10 Person | KS-MCC15C-03 | 740×575×1000(1070) | 15 Person | KS-MCC20C-03 | 740×575×1000(1070) | 20 Person |

- Storage Box (6 person) 2ea
- Lock (Key)

- Storage Box (6 person) 2ea
- Lock (Key)

- Storage Box (6 person) 2ea
- Lock (Key)

수술실용
Operating Room
Appliances간호차량
Nursing
Appliances진료실용
Consulting Room
Appliances공급실용
Consultation Room
Appliances의료용 침대
Hospital Bed카트
Carts모니터 카트
Monitor Carts수액 장아대
IV Pole Stand모니터 암
Monitor Arm물 이송도
Water Mount물 이송도
Pole Mount악세서리
AccessoryABS 카트
ABS Cart주문형 카트
Order Type Carts



MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE)

| Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark |
|----------|--------------------|----------|-----------|--------------------|-----------|-----------|--------------------|-----------|
| KS-MCC8C | 610×575×1000(1070) | 8 Person | KS-MCC12C | 610×575×1000(1070) | 12 Person | KS-MCC16C | 610×575×1000(1070) | 16 Person |

- Standard
- Lock (Key)
- Cassette Type

- Standard
- Lock (Key)
- Cassette Type

- Standard
- Lock (Key)
- Cassette Type



MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE)

| Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|--------|--------------|--------------------|--------|--------------|--------------------|--------|
| KS-MCC8C-01 | 610×575×1000(1070) | — | KS-MCC12C-01 | 610×575×1000(1070) | — | KS-MCC16C-01 | 610×575×1000(1070) | — |

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)
- Cassette Type

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)
- Cassette Type

- Utility Container
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)
- Cassette Type



MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE)

| Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|----------|--------------|--------------------|-----------|--------------|--------------------|-----------|
| KS-MCC8C-02 | 610×575×1000(1070) | 8 Person | KS-MCC12C-02 | 610×575×1000(1070) | 12 Person | KS-MCC16C-02 | 610×575×1000(1070) | 16 Person |

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)
- Cassette Type

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)
- Cassette Type

- Storage Box (6 person)
- Waste Holder
- Waste Box
- Needle Holder
- Lock (Key)
- Cassette Type



MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE) MEDICATION CART (610,CASSETTE)

| Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark | Model | Size (L×D×H / mm) | Remark |
|-------------|--------------------|----------|--------------|--------------------|-----------|--------------|--------------------|-----------|
| KS-MCC8C-03 | 610×575×1000(1070) | 8 Person | KS-MCC12C-03 | 610×575×1000(1070) | 12 Person | KS-MCC16C-03 | 610×575×1000(1070) | 16 Person |

- Storage Box (6 person) 2ea
- Lock (Key)
- Cassette Type

- Storage Box (6 person) 2ea
- Lock (Key)
- Cassette Type

- Storage Box (6 person) 2ea
- Lock (Key)
- Cassette Type



MEDICATION CART

| Model | Size (L×D×H / mm) | Remark |
|---------------|-------------------|--------|
| KS-SMCC25C-01 | 740×575×1000 | - |



MEDICATION CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|-----------|
| KS-MC2 | 630×520×950 | SUS(40인용) |

- 20/30인용 주문 제작 가능합니다.
- CA502~20EA



CASSETTE DELIVERY CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-CDLC-4B | 740×740×1400 | 4 BOX |



CASSETTE

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|-----------|
| KS-CB-10 | 563×330×250 | 10 Person |

※ Aluminium, A,B,S Box



CASSETTE

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|-----------|
| KS-CB-15 | 563×330×345 | 15 Person |

※ Aluminium, A,B,S Box



CASSETTE

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|-----------|
| KS-CB-20 | 563×330×440 | 20 Person |

※ Aluminium, A,B,S Box



CASSETTE

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|----------|
| KS-CB-8 | 463×330×250 | 8 Person |

※ Aluminium, A,B,S Box



CASSETTE

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|-----------|
| KS-CB-12 | 463×330×345 | 12 Person |

※ Aluminium, A,B,S Box



CASSETTE

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|-----------|
| KS-CB-16 | 463×330×440 | 16 Person |

※ Aluminium, A,B,S Box



BIN BOX

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|---------------|
| KS-AC-MB01 | 108×248×86 | A,B,S Plastic |

※ 칸 조절을 통하여 공간관리가 가능합니다.



BIN BOX

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|---------------|
| KS-AC-MB02 | 104×300×78 | A,B,S Plastic |

※ 칸 조절을 통하여 공간관리가 가능합니다.



ANESTHESIA CART

| Model | Size (L×D×H / mm) | Remark |
|--------------|--------------------|--------|
| KS-ATC-01-01 | 740×575×1070(1670) | — |



ANESTHESIA CART

| Model | Size (L×D×H / mm) | Remark |
|--------------|--------------------|--------|
| KS-ATC-01-02 | 740×575×1070(1670) | — |



ANESTHESIA CART

| Model | Size (L×D×H / mm) | Remark |
|--------------|--------------------|----------|
| KS-ATC-01-03 | 740×575×1070(1670) | 10 Hooks |



ICU CART

| Model | Size (L×D×H / mm) | Remark |
|------------|--------------------|--------|
| KS-ICUC-01 | 740×575×1000(1070) | — |



ANESTHESIA CART

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-ANC-01 | 640×530×1000~1485 | — |



KS-TRTC-5D



KS-TRTC-5D-01

SUPER PLUS CART

| Model | Size (L×D×H / mm) | Remark |
|---------------|--------------------|-----------|
| KS-TRTC-5D | 780×530×1000(1070) | Door Type |
| KS-TRTC-5D-01 | 780×530×1000(1070) | Open Type |



UTILITY CART

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------------|--------|
| KS-ANC-01 | 640×530×1000(1070~1485) | — |

FEATURE & USAGE

1. 장기간 입원 및 불편한 환자의 침대 곁에서 모발의 세척 및 건조가 가능하도록 한 보조 장치.
2. 내부 세척시 가능한 물탱크엔 최대 36Liter를 채워 2~3명에게 사용할 수 있음.
3. 탱크 내에서 물의 일정한 온도를 유지 할 수 있는 Digital 온도 조절장치 부착 및 전기 2차 안전 S/W 부착.
4. 절수기능이 뛰어난 샤워기 및 세척 후 바로 건조 할 수 있는 Hair Dryer 내장 되어 있음.
5. 사용한 오수를 처리하기 위한 저장 탱크가 하부에 부착되어 있음.
6. ATOM(Japan)사에서 직 수입된 Basin의 회전 및 높이 조절이 가능.

SPECIFICATION

- All Stainless Polish steel 304 Construction
- Power : 220V, 60Hz Code 4M
- Control Panel : Digital Panel, Main S/W, Heating S/W, Pump S/W
- Motor Pump : 1/8 Hp, Water Proof.
- Higher quality Brake Caster 4EA

ACCESSORIES



▷ Wash Basin(Atom)

▷ Hair Dryer

HAIR SHAMPOOING MACHINE

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-HSM1 | 770×470×850 | — |



WHEELCHAIR

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-KCA100 | 540×700×1040(875) | 병원용/스틸 |



영양이란 올랐을때

WALKER

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|------------|
| KS-WK01 | 750×450×950 | STANDARD |
| KS-O155H | 750×450×950 | TWO BUTTON |

◆ SPECIFICATION

- Stainless Steel Pipe 304 Construction
- Hi-Low Hand Control (3 Point)
- 영양이 받은 사용하지 않을 경우에는 뒤로 세워 놓을 수 있음
- 4" Swivel Caster 4EA (Brake 2EA)
- Accessories : LV Pole.

HAIR SHAMPOOING BASKET

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|---------|
| KS-HS-PA | 450×440×85 | 물조리개 포함 |

◆ SPECIFICATION

- Hair Shampooing Machine을 대체 가능
- 비교적 간단하게 머리를 감겨줄 수 있음
- 누구나 쉽게 사용할 수 있음
- 경제적인 가격
- 선물용 은 기증용으로도 좋음



R-ABS

UTILITY CART (DRAWER TYPE)

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|----------|
| KS-UC2S-2D | 665×500×950 | 2 Drawer |



Basket Drawer



Dust Box

ACCESSORIES

- Basket : 560×410×55(KS-AB-1)
- Basket : 560×410×125(KS-AB-2)
- Option : Dust Box

SPECIFICATION

- All Stainless Steel 304 Construction.
- Ball bearing drawer rail.
- ABS Basket in drawer.
- Length and breadth partition(Each a unit)
- Higher quality Caster 4EA
- Inner shelf size : 615×445×10



R-ABS

UTILITY CART (DRAWER TYPE)

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|----------|
| KS-UC3S-1D | 665×(500+350)×950 | 1 Drawer |



Basket Drawer



Sliding Plate

ACCESSORIES

- ABS Basket : 560×410×125(KS-AB-2)

SPECIFICATION

- All Stainless Steel 304 Construction.
- Ball bearing drawer rail.
- ABS Basket in drawer.
- Length and breadth partition(Each a unit)
- Higher quality Caster 4EA
- Inner shelf size : 615×445×10



UTILITY CART (A,B,S DRAWER)

| Model | Size (L×D×H / mm) |
|--------------|-------------------|
| KS-AUC2S-1DS | 565×500×950 |
| KS-AUC2S-1DL | 665×500×950 |

* R-ABS 적용



UTILITY CART (A,B,S DRAWER)

| Model | Size (L×D×H / mm) |
|--------------|-------------------|
| KS-AUC3S-1DS | 565×500×950 |
| KS-AUC3S-1DL | 665×500×950 |

* R-ABS 적용



UTILITY CART (A,B,S DRAWER)

| Model | Size (L×D×H / mm) |
|--------------|-------------------|
| KS-AUC2S-2DS | 565×500×950 |
| KS-AUC2S-2DL | 665×500×950 |

* R-ABS 적용



UTILITY CART (A,B,S DRAWER)

| Model | Size (L×D×H / mm) |
|--------------|-------------------|
| KS-AUC2S-H2D | 755×500×950 |

* R-ABS 적용



UTILITY CART (A,B,S DRAWER)

| Model | Size (L×D×H / mm) |
|--------------|-------------------|
| KS-AUC3S-H2D | 750×500×950 |

* R-ABS 적용



UTILITY CART (A,B,S DRAWER)

| Model | Size (L×D×H / mm) |
|--------------|-------------------|
| KS-AUC2S-H4D | 750×500×950 |

* R-ABS 적용



UTILITY CART

| Model | Size (L×D×H / mm) |
|-----------|-------------------|
| KS-UC2S-S | 650×450×900 |
| KS-UC2S-M | 750×450×900 |
| KS-UC2S-L | 950×450×900 |

* All Stainless Steel 304 Construction
 * 4" Swivel Caster 4EA (Brake 2EA)



UTILITY CART

| Model | Size (L×D×H / mm) |
|-----------|-------------------|
| KS-UC3S-S | 650×450×900 |
| KS-UC3S-M | 750×450×900 |
| KS-UC3S-L | 950×450×900 |

* All Stainless Steel 304 Construction
 * 4" Swivel Caster 4EA (Brake 2EA)



DRESSING CART

| Model | Size (L×D×H / mm) |
|-----------|-------------------|
| KS-DC2S-S | 650×450×900 |
| KS-DC2S-M | 750×450×900 |
| KS-DC2S-L | 950×450×900 |

* All Stainless Steel 304 Construction
 * 4" Swivel Caster 4EA (Brake 2EA)



DRESSING CART

| Model | Size (L×D×H / mm) |
|-----------|-------------------|
| KS-DC3S-S | 650×450×900 |
| KS-DC3S-M | 750×450×900 |
| KS-DC3S-L | 950×450×900 |

* All Stainless Steel 304 Construction
 * 4" Swivel Caster 4EA (Brake 2EA)



DRESSING CART (BANDING TYPE)

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-DC2SH-S | 650×450×950 |
| KS-DC2SH-M | 750×450×950 |
| KS-DC2SH-L | 950×450×950 |

* All Stainless Steel 304 Construction
 * One Body Pipe Frame Construction
 * 4" Swivel Caster 4EA (Brake 2EA)
 * R-ABS 적용



DRESSING CART (BANDING TYPE)

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-DC3SH-S | 650×450×950 |
| KS-DC3SH-M | 750×450×950 |
| KS-DC3SH-L | 950×450×950 |

* All Stainless Steel 304 Construction
 * One Body Pipe Frame Construction
 * 4" Swivel Caster 4EA (Brake 2EA)
 * R-ABS 적용



TREATMENT CART

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-TC2S1-S | 550×450×900 |
| KS-TC2S1-M | 650×450×900 |

- All Stainless Steel 304 Construction
- 4" Swivel Caster 4EA (Brake 2EA)



TREATMENT CART

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-TC2S2-S | 550×450×900 |
| KS-TC2S2-M | 650×450×900 |

- All Stainless Steel 304 Construction
- 4" Swivel Caster 4EA (Brake 2EA)
- 상기 모델은 원하는 규격으로 주문 생산이 가능합니다.



TREATMENT CART

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-TC2S2-S | 750×450×900 |
| KS-TC2S2-M | 850×450×900 |
| KS-TC2S2-L | 950×450×900 |

- All Stainless Steel 304 Construction
- 4" Swivel Caster 4EA (Brake 2EA)



TREATMENT CART

| Model | Size (L×D×H / mm) |
|------------|-------------------|
| KS-TC3S2-S | 750×450×900 |
| KS-TC3S2-M | 850×450×900 |
| KS-TC3S2-L | 950×450×900 |

- All Stainless Steel 304 Construction
- 4" Swivel Caster 4EA (Brake 2EA)



TREATMENT CART (BANDING TYPE)

| Model | Size (L×D×H / mm) |
|-------------|-------------------|
| KS-TC2S2H-S | 750×450×950 |
| KS-TC2S2H-M | 850×450×950 |
| KS-TC2S2H-L | 950×450×950 |

- All Stainless Steel 304 Construction
- 4" Swivel Caster 4EA (Brake 2EA)
- R-ABS 적용



TREATMENT CART (BANDING TYPE)

| Model | Size (L×D×H / mm) |
|-------------|-------------------|
| KS-TC3S2H-S | 750×450×950 |
| KS-TC3S2H-M | 850×450×950 |
| KS-TC3S2H-L | 950×450×950 |

- All Stainless Steel 304 Construction
- 4" Swivel Caster 4EA (Brake 2EA)
- R-ABS 적용




SIDE TABLE

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|------------------------|
| KS-ST-01 | 800×350×775 | 2" Caster, Wood, Steel |
| KS-ST-02 | 800×350×775 | 1" Caster, H45 |



SIDE TABLE (C-ARM)

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|----------|
| KS-ST-03 | 700×400×775 | BAKE 15T |



URINE CONTAINER RACK (HOLDER TYPE)

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|----------|
| KS-UCR01 | 100×220×260 | 1 Holder |
| KS-UCR02 | 200×220×260 | 2 Holder |
| KS-UCR03 | 300×220×260 | 3 Holder |

- All Stainless Steel 304 Construction
- Pipe Frame & Joint Type



FECES PEN RACK (HOLDER TYPE)

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|----------|
| KS-FPR01 | 100×220×330 | 1 Holder |
| KS-FPR02 | 200×220×330 | 2 Holder |
| KS-FPR03 | 300×220×330 | 3 Holder |

- All Stainless Steel 304 Construction
- Pipe Frame & Joint Type



SIDE TABLE (ARM SUPPORT)

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|-----------------|
| KS-ST-04 | 460×570×775 | 매트리스 600×150×40 |



BED PEN RACK

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|----------|
| KS-BPR03 | 500×420×800 | 3 Holder |
| KS-BPR05 | 850×420×800 | 5 Holder |
| KS-BPR08 | 1100×420×800 | 8 Holder |

- All Stainless Steel 304 Construction
- Drain Hole & Tray
- Wall Mounted Type



LINEN HAMPER CART

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-LHCF | 600×600×850 | — |
| KS-LHCF-1 | 900×600×850 | — |
| KS-LHCF-2 | 1200×600×850 | — |

- Stainless Steel Pipe 304 Frame
- 4" Swivel Caster 4EA
- Accessories : Hamper bag 1EA
- Hamper bag 의 색상은 자체 수급에 따라 바뀔 수 있음



LINEN HAMPER CART

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-LHCR | φ600×850 | — |

- Stainless Steel Pipe 304 Frame
- 4" Swivel Caster 3EA
- Accessories : Hamper bag 1EA
- Hamper bag 의 색상은 자체 수급에 따라 바뀔 수 있음



LINEN HAMPER CART

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|------------|
| KS-LHCF-C | 500×500×850 | Cover Type |

- Stainless Steel Pipe 304 Frame
- Cover : Acrylic 5.0mm & Handle
- 3" Swivel Caster 4EA
- Accessories : Hamper bag 1EA
- Hamper bag 의 색상은 자체 수급에 따라 바뀔 수 있음
- Pedal Type

LINEN HAMPER BAG (햄퍼포)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-LHB-01 | 500×500 | — |
| KS-LHB-02 | 600×600 | — |

- OPTION-방수포(오렌지, 녹색)



LINEN HAMPER CART

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-LHCD | 1000×450×850 | Double |

- Stainless Steel Pipe 304 Frame
- Stainless Steel Shelf Plate
- Accessories : Hamper bag 2EA
- Hamper bag 의 색상은 자체 수급에 따라 바뀔 수 있음



STEP FOOT STOOL

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|----------|
| KS-FS1-1 | 430×300×120 | One Step |
| KS-FS1-2 | 430×300×180 | One Step |
| KS-FS1-3 | 430×300×250 | One Step |

- All Stainless Steel 304 Construction
- Non Slip rubber top surface
- Non Slip rubber Parking



STEP FOOT STOOL

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|----------|
| KS-FS2-1 | 430×500×500 | One Step |
| KS-FS2-2 | 700×500×500 | One Step |

- All Stainless Steel 304 Construction
- Non Slip rubber top surface
- Non Slip rubber Parking



STEP FOOT STOOL

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-FS-3 | 430×500×470 | — |

- All Stainless Steel 304 Construction



CHART HOLDER CART

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-CHC530 | 530×335×1000 | 40매 |
| KS-CHC630 | 630×335×1000 | 60매 |
| KS-CHC830 | 830×335×1000 | 80매 |

- All Stainless Steel 304 Construction
- 4" Swivel Caster 4EA



CHART HOLDER CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-CHC530S | 530×335×1200 | 40매 |
| KS-CHC630S | 630×335×1200 | 60매 |
| KS-CHC830S | 830×335×1200 | 80매 |

- All Stainless Steel 304 Construction
- 4" Swivel Caster 4EA



BASKET CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-BC1 | 600×450×900 | — |



BASKET CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-BC2 | 600×450×900 | — |



BASKET CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-BC3 | 700×500×950 | — |



OXYGEN TANK CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-OTC10-A | — | φ120 |
| KS-OTC10-B | — | φ160 |

* All SUS



OXYGEN TANK HOLDER

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-OTH120 | φ120×440 | — |



OXYGEN TANK HOLDER

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-OTH160 | φ160×370 | — |



OXYGEN TANK CART

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-OTC20 | For 20kg | — |

* SUS, Steel



OXYGEN TANK CART

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-OTC40 | For 40kg | — |

* SUS, Steel



SCREEN

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|---------|
| KS-1PS | 900×440×1800 | 1 Panel |



SCREEN

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|---------|
| KS-2PS | 1800×400×1800 | 2 Panel |



SCREEN

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|---------|
| KS-3PS | 2700×400×1800 | 3 Panel |

- ※ All Stainless Steel 304 Construction
- ※ 2" Swivel Caster
- ※ Screen Cloth : Pure cotton



환자운반카 (유압식)

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|----------|
| KS-HDSC200 | 1825×615×540~870 | 무게 55 kg |

- ※ 승하강 유압방식
- ※ 통각도 조절방식



환자운반카 (TOTAL LOCK)

| Model | Size (L×D×H / mm) | Remark |
|---------------|-------------------|--------|
| KS-HDSC-TW400 | 955×700×450~895 | - |

- ※ 승하강 유압방식
- ※ 통각도 조절방식

수술실용
Operating Room
Appliances간호실용
Nursing
Appliances진료실용
Consulting Room
Appliances환자실용
Supply Room
Appliances의료용 침대
Hospital Bed

카트

모니터 카트

수액 관리대

IV Pole Stand

모니터 암

Wall Mount

Pole Mount

Accessories


ABS Cart

수문용 카트

Order Type Carts



1 걸이대



2 카테터걸이



3 잠금장치



위 내시경장 (ENDOSCOPIC CABINET)

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|----------|
| KS-ESC-2 | 700×450×2100 | 2 Holder |
| KS-ESC-4 | 900×450×2100 | 4 Holder |
| KS-ESC-6 | 1200×450×2100 | 6 Holder |

대장 내시경장 (ENDOSCOPIC CABINET)

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|----------|
| KS-CSC-2 | 700×450×2300 | 2 Holder |
| KS-CSC-4 | 900×450×2300 | 4 Holder |
| KS-CSC-6 | 1200×450×2300 | 6 Holder |

특징 및 구성

1. 부식성에 뛰어난 스테인레스 304 재질로 제작됩니다.
2. 내시경 보관의 안전성과 효율성을 유지하기 위하여 별도의 보관을 위하여 제작됩니다.
3. 세척 후 건조와 보관을 할 수 있도록 합니다. Endoscope를 안전하고 편리하게 좁은 공간에서 더 많이 보관하도록 제작됩니다.
4. 별도의 Fan을 부착하여 내부의 습도를 조절합니다.
5. 5T 안전강화유리 사용으로 내부의 상태를 확인합니다.
6. 미닫이문에 달린 잠금 장치로 분실에 대비합니다.

위 내시경 회전장
ENDOSCOPIC CABINET(ROTATE TYPE)

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|-----------|
| KS-ESC-R6 | 700×700×2100 | 6 Holder |
| KS-ESC-R8 | 700×700×2100 | 8 Holder |
| KS-ESC-R10 | 700×700×2100 | 10 Holder |

대장 내시경 회전장
ENDOSCOPIC CABINET(ROTATE TYPE)

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|-----------|
| KS-CSC-R6 | 700×700×2300 | 6 Holder |
| KS-CSC-R8 | 700×700×2300 | 8 Holder |
| KS-CSC-R10 | 700×700×2300 | 10 Holder |

특징 및 구성

1. 부식성에 뛰어난 스테인레스 304 재질로 제작됩니다.
2. 내시경 보관의 안전성과 효율성을 유지하기 위하여 별도의 보관을 위하여 제작됩니다.
3. 세척 후 건조와 보관을 할 수 있도록 합니다. Endoscope를 안전하고 편리하게 좁은 공간에서 더 많이 보관하도록 제작됩니다.
4. Holder는 회전할 수 있도록 되어 있습니다.
5. 별도의 Fan을 부착하여 내부의 습도를 조절합니다.
6. 5T 안전강화유리 사용으로 내부의 상태를 확인합니다.
7. 미닫이문에 달린 손잡이에 잠금 장치를 적용합니다.



1 걸이대



2 홀더



3 홀더





ACCESSORIES STAND

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-AS-02 | Φ500×1500~2200 | - |

※ 2" Caster, SUS 304



ACCESSORIES STAND

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-AS-03 | Φ500×1500×2200 | - |

※ 2" Caster, SUS 304



ENDOSCOPIC DRYING RACK (WOOD)

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-EH-W1 | 75×93×103 | U Type |

※ Plastic
※ 단독 부착 시 조임볼트 제외



ENDOSCOPIC DRYING RACK (WOOD)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-NEH-W1 | 120×40×40 | R Type |

※ Plastic (Holder) + Steel (Bracket)
※ R Type Holder 의 단독부착 할 수 없음



ENDOSCOPIC DRYING RACK STAND

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-EH-4 | Φ550×1200~1800 | 4구 |

※ 3" Caster ABS Base, Stainless, Aluminium, Rubber
※ 높낮이 조절가능



ENDOSCOPIC DRYING RACK STAND

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-EH-6 | Φ630×1200~1900 | 6구 |

KS-EH-8 Φ630×1200~1900 8구
KS-EH-10 Φ630×1200~1900 10구
※ 3" Caster ABS Base, Stainless, Aluminium, Rubber
※ 높낮이 조절가능



ENDOSCOPIC DRYING RACK (WOOD)

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|----------|
| KS-EH-W2 | W=380 | 2 Holder |

※ 간격은 기본사이즈입니다.
※ 요구에 따라 사이즈 제작 가능합니다.
※ Bar : Aluminium



ENDOSCOPIC DRYING RACK (WOOD)

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|----------|
| KS-EH-W3 | W=710 | 3 Holder |


※ 간격은 기본사이즈입니다.
※ 요구에 따라 사이즈 제작 가능합니다.
※ Bar : Aluminium



ENPOSCOPIC LENS SAFER

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-SH-1 | Φ25~20 | - |

- Endoscope Lens Safer
- Rubber (Holder)
- U or R Type과 함께 사용



ENPOSCOPIC LENS SAFER

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-SH-S1 | W=100 | - |


- Endoscope Lens Safer Bracket/Holder
- Rubber (Holder) + Steel (Bracket)
- U or R Type과 함께 사용



ENPOSCOPIC LENS SAFER

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-SH-S2 | W=245 | - |

- Endoscope Lens Safer Bracket 2 Holder
- Rubber (Holder) + Steel (Bracket)
- U or R Type과 함께 사용



FILM CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-FC1 | 650×500×850 | SUS |

- 5" Caster, Corner Safety Rubber



FILM CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-FC2 | 750×450×950 | SUS |


- 5" Caster



APRON STORAGE RACK

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|-----------|
| KS-ASR | 500×500×1450 | 10 Holder |

- 3" Caster, Stainless Steel 304
- Φ10 Bar & Steel Frame



X-RAY ILLUMINATOR

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-XR100 | 485×500×39 | 고급형 |
| KS-XR200 | 850×500×75 | 보급형 |
| KS-XR300 | 1210×570×40 | 고급형 |
| KS-XR400 | 1590×500×40 | 고급형 |



FILM CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-FC3 | 1200×500×1200 | SUS |

- 2.5" Caster, Corner Safety Rubber



FILM CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-FC4 | 700×550×1150 | Steel |

- 5" Caster, Rubber



OPERATING STOOL

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-271050 | 370×460×880~960 | 진찰용 A형 |



OPERATING STOOL

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-263330 | φ370×510~700 | 진찰용 B형 |



OPERATING STOOL

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-262540 | 380×440×770~950 | 진찰용 D형 |

* 검사발



OPERATING STOOL

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-270050 | 380×440×770~950 | 진찰용 D형 |

* 오리발



OPERATING STOOL

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|-------------|
| KS-264820 | φ370×510~760 | 진찰용 C형, 이물식 |
| KS-264820C | φ370×510~760 | 진찰용 C형, 고정식 |



OPERATING STOOL

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-267850 | 370×545×780~900 | 진찰용 O형 |



OPERATING STOOL

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-DR-400 | 610×550×900~1170 | 제도용 |



OPERATING STOOL

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-312230 | φ370×510~700 | 소아용 |



WIRE SHELF CART

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-WC2H-S | 910×460×950 | — |
| KS-WC2H-L | 910×610×950 | — |

- Stainless Steel 304 Pipe Frame
- 5" Swivel Caster 4EA
- Adjustable Wire Shelves (30mm)



WIRE SHELF CART

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-WC3-S | 910×460×950 | — |
| KS-WC3-L | 910×610×950 | — |

- Stainless Steel 304 Pipe Frame
- 5" Swivel Caster 4EA
- Adjustable Wire Shelves (30mm)



WIRE SHELF CART

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-WC2-S | 910×460×950 | — |
| KS-WC2-L | 910×610×950 | — |

- Stainless Steel 304 Pipe Frame
- 5" Swivel Caster 4EA
- Adjustable Wire Shelves (30mm)



WIRE SHELF CART

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-WC3H-S | 910×460×950 | — |
| KS-WC3H-L | 910×610×950 | — |

- Stainless Steel 304 Pipe Frame
- 5" Swivel Caster 4EA
- Adjustable Wire Shelves (30mm)



WIRE SHELF

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-W910S | 910×460×1800 | — |
| KS-W910L | 910×610×1800 | — |
| KS-W1220S | 1220×460×1800 | — |
| KS-W1220L | 1220×610×1800 | — |
| KS-W1520S | 1520×460×1800 | — |
| KS-W1520L | 1520×610×1800 | — |
| KS-W1820S | 1820×460×1800 | — |
| KS-W1820L | 1820×610×1800 | — |

- Stainless Steel 304 Pipe Frame
- Height adjustable all shelf
- Options : Guardrail



WIRE MASH (OPTION)

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-MASH-01 | 추후협의 | — |



WIRE SHELF CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-WC910S | 910×460×1800 | — |
| KS-WC910L | 910×610×1800 | — |
| KS-WC1220S | 1220×460×1800 | — |
| KS-WC1220L | 1220×610×1800 | — |
| KS-WC1520S | 1520×460×1800 | — |
| KS-WC1520L | 1520×610×1800 | — |
| KS-WC1820S | 1820×460×1800 | — |
| KS-WC1820L | 1820×610×1800 | — |

- Stainless Steel 304 Pipe Frame
- Height adjustable all shelf
- Caster : 5" 4EA
- Options : Guardrail



PLASTIC DOUBLE CART

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|---------|
| KS-PDSC | 750×460×850 | 하중 60kg |

- Double plate available
- Plastic body is light in weight and gives easy moving
- Removable handle



PLASTIC TRIPLE CART

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|---------|
| KS-PTSC | 750×460×850 | 하중 60kg |

- Triple plate available
- Plastic body is light in weight and gives easy moving
- Removable handle



01
사이드프레임
기둥과 기둥사이를 고정
시켜주는 악세사리



02
난간
선반에 보관중인 제품의
낙하를 방지



03
디바이더
분리용 칸막이, 제품을 분리
및 보관에 용이



04
손잡이 핸들
이동의 편리성을 위해
손잡이 결합가능



05
S 고리
선반과 선반을 직각 또는
일렬로 배열할 때 사용



06
평평형브라켓
선반을 일렬로 배열시 기둥과
기둥을 결합 시킬 때



07
높이조정볼트
메탈랙의 높낮이를 조절할
수 있는 조정 나사




08
슬리브
기둥에 고정되어 선반과
결합시킴

SLIDING METAL RACK SYSTEM

| Model | Size | | | Remark | |
|---------|------|-----|------|--|------|
| | W | D | H | | |
| KS-SMRS | 610 | 460 | 1500 | 단수/높이/ 레일길이는 공간에 맞추어 선택 조절 가능 | |
| | 760 | | | | |
| | 910 | | | | |
| | 1070 | 610 | | | 1800 |
| | 1220 | | | | 1900 |
| | 1520 | | 2200 | | |
| | 1820 | 450 | | | |
| | 750 | | | | |
| | 900 | | | | |
| | 1200 | | | | |

= Chromium, Stainless, Aluminium



우레탄바퀴
메탈랙/카트의 이동을 간편하게 할
수 있도록 바퀴를 장착할 수 있습니다.
무거운 물건을 옮기고 자주 이동 할
경우에는 5인치 바퀴를 사용 하시는
것이 좋습니다.



TRUCK CART

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-TC1 | 900×650×850 | - |
| KS-TC1-1 | 1200×650×850 | - |

※ Stainless Steel 304 Construction
※ 5" Swivel Caster 4EA



TRUCK CART

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-TC2 | 900×650×850 | - |
| KS-TC2-1 | 1200×650×850 | - |

※ Stainless Steel 304 Construction
※ 5" Swivel Caster 4EA



DANPRA(OPTION)

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-DP-01 | 추후협의 | - |



WORK TABLE

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-BWT01 | 2300×1100×900 | - |
| KS-BWT02 | 2000×1100×900 | - |

※ Body : Stainless Steel Plate 304,
(Top Plate Support : 12.0T Wood Plate)
※ Higher Quality Caster 4EA,
※ Basket Composition,
※ A.B.S Basket
: 560X410X55H & 560X410X125H
※ SUS Wire Basket
: 560X410X125H & 560X410X190H



LINEN SERVICE CART

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-LSC1-M | 1200×600×1600 | - |
| KS-LSC1-L | 1500×600×1600 | - |

※ Stainless Steel 304 Construction
※ 5" Swivel Caster 4EA



LINEN SERVICE CART

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-LSC2-M | 1200×800×1200 | - |
| KS-LSC2-L | 1500×800×1200 | - |

※ Stainless Steel 304 Construction
※ 5" Swivel Caster 4EA



DELIVERY CART

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-DVC1 | 1200×620×1030 | - |
| KS-DVC2 | 1400×620×1030 | - |

※ Stainless Steel 304 Construction
※ 5" Swivel Caster 4EA



WORK TABLE

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-WRT01 | 1500×700×900 | - |

※ Stainless Steel Plate 304 Construction
※ Higher quality Caster 4EA



WORK TABLE

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-WRT02 | 1500×700×900 | - |

※ Stainless Steel Plate 304 Construction
※ Higher quality Caster 4EA



SEALING TABLE

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-SST01 | 900×600×1800 | - |

- Shelf & Drawer : 304 Stainless Steel Plate
- Support : AL Profile
- SUS Wire Basket (400×300×150H)
- Q'TY : 3EA
- Castor : High Quality 5 inch 4EA

SEALING STAND

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-SS01 | 800×380×700 | - |

- Body : 304 Stainless Steel Plate
- Support : 304 Stainless Steel Pipe
- Seal Guard : Stainless Bar Φ8 (2EA)
- Seal Roller : MC Plastic



SAFETY MEDICINE CABINET

| Model | Size (L×D×H / mm) | Remark |
|-----------|----------------------------------|--------|
| KS-C119EN | 외부 412×363×360 내부 320×230×260 | 디지털 락 |

- 용적 : 19ℓ
- 중량 : 33kg



SAFETY MEDICINE CABINET

| Model | Size (L×D×H / mm) | Remark |
|-----------|----------------------------------|-----------|
| KS-C119EK | 외부 412×363×360 내부 320×230×260 | 디지털 락+키 락 |

- 용적 : 19ℓ
- 중량 : 33kg



내부구성 및 특징

- 별도 잠금장치가 부착된 속문
- 튼튼한 후레임 구조로 연기, 열기 차단
- 좌우, 2방향 잠금핀으로 방도성 강화
- 문안 쪽 키행거 4개
- 특수 표면처리 : 긁힘과 녹방지



EN (디지털락)



EK (디지털락 + 키락)

디지털락 EDL88 특징

- FND창 : 숫자 및 기능 표시
- 비밀번호 : 4~16자리까지 설정
- 이중잠금기능 : 2개의 비밀번호를 입력해야 문이 열림
- 허수기능 : 비밀번호 앞뒤에 허수를 넣어 사용자 비밀번호 노출 방지
- 마비기능 : 6회이상 틀린번호 입력시 마비
- 감충기능 : 비밀번호와 소리 숨김
- 경보기능 : 금고가 45도 이상 기울거나 충격이 가해지면 경보를 발생
- 건전지 : 1.5V알카라인 AA건전지 4개
- 건전지교체 알림 : FND창에 LOW-BAT 표시



DUST BOX CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-DBC | Φ500×800 | - |

- Stainless Steel 304 Construction
- 4" Swivel Caster 4EA



LINEN SERVICE CART

| Model | Size (L×D×H / mm) | Remark |
|-----------|--------------------|--------|
| KS-LSC3-1 | 740×575×1000(1070) | - |



FOOD CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-----------------------------|--------|
| KS-FD1 | 식판 사이즈에 따라 변경 40/48 Persons | - |

- Stainless Steel 304 Construction
- 5" Swivel Caster 4EA




SAFETY MEDICINE CABINET

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-SMC-01 | 1206×654×1799 | 3Panel |

- 내부: 1047×487×1536 (W×D×H)
- 선반 3개
- 서랍 2개
- 중량 670KG
- 잠금방식 : 디지털락, 키(Key)





MULTI BIN (PLASTIC)


| Model | Size (L×D×H / mm) | Remark |
|-------------------|-------------------|--------|
| KS-CA500 | 82×162(150)×41 | — |
| ※ 서랍규격: 73×146×32 | | |
| ※ 칸막이규격: 69×26 | | |

MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|-------------------|-------------------|--------|
| KS-CA501 | 159×162(150)×41 | — |
| ※ 서랍규격: 73×146×32 | | |
| ※ 칸막이규격: 69×26 | | |

MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|--------------------|-------------------|--------|
| KS-CA501-2 | 216×207(189)×63 | — |
| ※ 서랍규격: 102×185×53 | | |
| ※ 칸막이규격: 94×38 | | |



MULTI BIN (PLASTIC)


| Model | Size (L×D×H / mm) | Remark |
|--------------------|-------------------|--------|
| KS-CA504-2 | 267×345(328)×99 | — |
| ※ 서랍규격: 126×322×45 | | |
| ※ 칸막이규격: 119×30 | | |

MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|--------------------------------|-------------------|--------|
| KS-CA505 | 267×345(328)×99 | — |
| ※ 서랍규격: 257×322×45 | | |
| ※ 칸막이규격: 가로 247×29 / 세로 315×29 | | |

MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|------------------------------|-------------------|--------|
| KS-CA506 | 267×345(328)×99 | — |
| ※ 서랍규격: 257×322×90 | | |
| ※ 칸막이규격: 가로 247×60 세로 315×60 | | |



MULTI BIN (PLASTIC)


| Model | Size (L×D×H / mm) | Remark |
|--------------------|-------------------|--------|
| KS-CA502 | 260×235(218)×72 | — |
| ※ 서랍규격: 122×212×62 | | |
| ※ 칸막이규격: 114×45 | | |

MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|--------------------|-------------------|--------|
| KS-CA503 | 260×235(218)×72 | — |
| ※ 서랍규격: 250×212×62 | | |
| ※ 칸막이규격: 없음 | | |

MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|--------------------|-------------------|--------|
| KS-CA504 | 267×345(328)×99 | — |
| ※ 서랍규격: 126×322×90 | | |
| ※ 칸막이규격: 119×75 | | |



MULTI BIN (PLASTIC)


| Model | Size (L×D×H / mm) | Remark |
|--------------------|-------------------|--------|
| KS-CA507 | 400×235(217)×295 | — |
| ※ 서랍규격: 122×212×62 | | |
| ※ 칸막이규격: 없음 | | |

MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|-------------------|-------------------|--------|
| KS-CA508 | 330×164×227 | — |
| ※ 서랍규격: 73×146×32 | | |
| ※ 칸막이규격: 없음 | | |

MULTI BIN (PLASTIC)


| Model | Size (L×D×H / mm) | Remark |
|--------------------|-------------------|--------|
| KS-CA509 | 344×205(191)×270 | — |
| ※ 서랍규격: 102×185×53 | | |
| ※ 칸막이규격: 94×38 | | |



MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|-------|-------------------|--------|
|-------|-------------------|--------|


※ 서랍규격: 124×326×102
※ 칸막이규격: 120×85



MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|-------|-------------------|--------|
|-------|-------------------|--------|


※ 서랍규격: 128×326×49
※ 칸막이규격: 120×35



MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|-------|-------------------|--------|
|-------|-------------------|--------|



※ 서랍규격: 258×326×49
※ 칸막이규격: 가로 250×35 / 세로 319×35



MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|-------|-------------------|--------|
|-------|-------------------|--------|



※ 서랍규격: 447×322×91
※ 칸막이규격: 가로 438×70 / 세로 314×70



MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|-------|-------------------|--------|
|-------|-------------------|--------|


※ 서랍규격: 447×322×60
※ 칸막이규격: 가로 438×45 / 세로 314×45



MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|-------|-------------------|--------|
|-------|-------------------|--------|


※ 서랍규격: 447×322×187
※ 칸막이규격: 가로 438×140 / 세로 314×140



MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|-------|-------------------|--------|
|-------|-------------------|--------|


※ 서랍규격: 258×326×102
※ 칸막이규격: 가로 250×84 / 세로 319×84



MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|-------|-------------------|--------|
|-------|-------------------|--------|



※ 서랍규격: 220×322×187
※ 칸막이규격: 210×140



MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|-------|-------------------|--------|
|-------|-------------------|--------|



※ 서랍규격: 220×322×91
※ 칸막이규격: 210×70



MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|-------|-------------------|--------|
|-------|-------------------|--------|


※ 서랍규격: 135×207×103
※ 칸막이규격: 128×84



MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|-------|-------------------|--------|
|-------|-------------------|--------|


※ 서랍규격: 155×322×144
※ 칸막이규격: 149×117




MULTI BIN (PLASTIC)

| Model | Size (L×D×H / mm) | Remark |
|-------|-------------------|--------|
|-------|-------------------|--------|

※ 서랍규격: 220×322×187
※ 칸막이규격: 210×140



- 세랍규격: CA505서랍과 동일
- 칸막이규격: 없음



| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-CA704 | 201×160×250 | — |



- * Aluminium, MDF, Plastic
- * CA704 : 3EA
- * CA703 : 8EA
- * CA702 : 18EA



- Aluminium, MDF, Plastic
- CA507 : 192EA



MULTI BIN CABINET (WOOD)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-CA1113 | 842×400×1600 | — |

- ※ CA514:14EA
- ※ CA514-2:7EA
- ※ CA515:7EA
- ※ CA516:14EA



MULTI BIN CABINET (WOOD)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-CA1114 | 1109×400×1600 | — |

- ※ CA514:2EA
- ※ CA514-2:2EA
- ※ CA515:8EA
- ※ CA516:8EA



MULTI BIN CABINET (WOOD)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-CA1101 | 1105×400×1100 | — |

- ※ CA504:10EA
- ※ CA504-2:10EA
- ※ CA505:10EA
- ※ CA506:10EA



MULTI BIN CABINET (WOOD)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-CA1115 | 954×400×1650 | — |

- ※ CA524:8EA
- ※ CA524-2:4EA
- ※ CA525:2EA
- ※ CA526:2EA



MULTI BIN CABINET (WOOD)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-CA1116 | 1414×400×1650 | — |

- ※ CA524:15EA



MULTI BIN CABINET (WOOD)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-CA1102 | 1375×400×1650 | — |

- ※ CA504:20EA
- ※ CA504-2:10EA
- ※ CA505:10EA
- ※ CA506:10EA



MULTI BIN CABINET (WOOD)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-CA1104 | 1420×460×1650 | — |

- ※ CA504:20EA
- ※ CA504-2:10EA
- ※ CA505:10EA
- ※ CA506:10EA



MULTI BIN CART

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-SS-33 | 600×520×1520 | — |

※ 인수에 따라 변형 제작가능
 ※ CA702:36EA
 ※ CA703:24EA
 ※ CA704:6EA



PASS BOX SYSTEM

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|-----------|
| KS-PBS36 | | 36 Person |

※ Sun 304
 ※ Numbr Key Unit
 ※ 간수 변형가능



TRAY RACK CART (SUS)

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-TRC | 1400×620×1320 | — |

※ 5" Caster, SUS 304, 60 Persons



U,V STERILIZER

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|---------|
| KS-US-050H | 500×410×510 | 250Watt |

※ CA524:15EA

MONITOR ARM



CHART CASE



BED MONITOR ARM

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|-------------|
| KS-BMA-01 | | MONITOR ARM |
| KS-BC-01 | | CHART CASE |

[3 MOTOR]



MOTORIZED PATIENT-CONTROL BED


| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-KY-001-1 | 2080×850×400 | 전동식 |

[2 MOTOR]




MOTORIZED PATIENT-CONTROL BED

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-KY-002-2 | 2080×850×400 | 전동식 |



3MOTOR GATCH BED

- All motorized, head & knee positions
- Detachable head and foot boards are made of polyester dressed polywood with steel reinforced plastic grips at both ends(Made-order production)
- Motor : AC 220v 50 ~ 60Hz DC 24V Dewert, Germany
- Caster : 4",Swevie Caster
- Finish : Electrostatic Powder Spray Coating



2MOTOR GATCH BED

- All motorized, head & knee positions
- Detachable head and foot boards are made of polyester dressed polywood with steel reinforced plastic grips at both ends(Made-order production)
- Motor : AC 220v 50 ~ 60Hz DC 24V Dewert, Germany
- Caster : 4",Swevie Caster
- Finish : Electrostatic Powder Spray Coating



[1 MOTOR]



MOTORIZED PATIENT-CONTROL BED

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-KY-002-3 | 2080×850×400 | 전동식 |



1MOTOR GATCH BED

- All motorized, head & knee positions
- Detachable head and foot boards are made of polyester dressed plywood with steel reinforced plastic grips at both ends(Made-to-order production)
- Motor : AC 220v 50 ~ 60Hz DC 24V Dewert, Germany
- Caster : 4", Swivel Caster
- Finish : Electrostatic Powder Spray Coating



[2 CRANK]



MOTORIZED PATIENT-CONTROL BED

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|------------|
| KS-KY-005 | 2080×850×500 | 2CRANK BED |



OPTIONAL ACCESSORY

- Dependent : 1EA
- Side rail : 2EA
- Mattress : 1EA
- Caster : 1SET(4"Swivel Caster 4EA)
- I.V.Pole : 1EA

SPECIFICATION

- Bed Top Frame : T 2.3mm Steel Plate
- Bed Top : Steel plate
- Head/Foot Board : P-P(Poly Propylene)
- Crank : Φ32 Steel Pipe, Universal joint crank
- HandleNylon resins 2EA (Slide Type)
- Bumper : T 25 X Φ100mm Nylon resins 4EA
- FinishElectrostatic Power Spray Coating (Use of the epoxy-polyester resins)




[1 CRANK]



HOSPITAL BED

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|------------|
| KS-KY-006 | 2080×850×500 | 1CRANK BED |



OPTIONAL ACCESSORY

- Dependent : 1EA
- Caster : 1SET(4*Swivel Caster 4EA)
- Side rail : 2EA
- I.V.Pole : 1EA
- Mattress : 1EA

SPECIFICATION

- Bed Top Frame : T 2.3mm Steel Plate
- Bed Top : Steel plate
- Head/Foot Board : P-P(Poly Propylene)
- Crank : $\Phi 32$ Steel Pipe, Universal joint crank
- Handle : Nylon resins 2EA (Slide Type)
- Bumper : T 25 X $\Phi 100$ mm Nylon resins 4EA
- Finish : Electrostatic Power Spray Coating (Use of the epoxy-polyester resins)




[1 CRANK]



소아용침대

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|------------|
| KS-KY-011 | 1880×850×500 | 1CRANK BED |



OPTIONAL ACCESSORY

- Dependent : 1EA
- Caster : 1SET(4*Swivel Caster 4EA)
- Side rail : 2EA
- I.V.Pole : 1EA
- Mattress : 1EA

SPECIFICATION

- Bed Top Frame : T 2.3mm Steel Plate
- Bed Top : Steel plate
- Head/Foot Board : P-P(Poly Propylene)
- Crank : $\Phi 32$ Steel Pipe, Universal joint crank
- Handle : Nylon resins 2EA (Slide Type)
- Bumper : T 25 X $\Phi 100$ mm Nylon resins 4EA
- Finish : Electrostatic Power Spray Coating (Use of the epoxy-polyester resins)




[1 CRANK]



소아용침대

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|------------|
| KS-KY-012 | 1880×850×500 | 1CRANK BED |



OPTIONAL ACCESSORY

- Dependent : 1EA
- Caster : 1SET(4'Swivel Caster 4EA)
- Side rail : 2EA
- LV.Pole : 1EA
- Mattress : 1EA

SPECIFICATION

- Bed Top Frame : T 2.3mm Steel Plate
- Bed Top : Steel plate
- Head/Foot Board : P-P(Poly Propylene)
- Crank : Φ32 Steel Pipe, Universal joint crank
- Handle : Nylon resins 2EA (Slide Type)
- Bumper : T 25 X Φ100mm Nylon resins 4EA
- Finish : Electrostatic Power Spray Coating (Use of the epoxy-polyester resins)



PILLOW

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-PL-01 | 추후 협의 | - |

= 사이즈 변경 가능

MATTRESS

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|-----------|
| KS-MAT01 | 1800×700×50 | 일반 (50T) |
| KS-MAT02 | 1870×880(830)×50 | 소아용 (50T) |

= 사이즈 변경 가능

ATTENDANT BED

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-KY-014 | 1600×550×250 | - |



EXAMINING TABLE

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-KY-025 | 1800×600(700)×650 | - |

SPECIFICATION

- Table Top : T11mm, Polywood with Sponge & Vinyl leather cover 40mm
- Leg : 30 X 50mm Steel Pipe



EXAMINING TABLE

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-KY-026 | 1800×600(700)×650 | - |

SPECIFICATION [1-HAND WORK]

- Table Top : T11mm, Polywood with Sponge & Vinyl leather cover 40mm
- Leg : 30 X 50mm Steel Pipe
- Side rail : 2EA



EXAMINING TABLE

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-KY-027 | 1800×600(700)×650 | - |

SPECIFICATION [1-HAND WORK]

- Table Top : T11mm, Polywood with Sponge & Vinyl leather cover 40mm
- Leg : 30 X 50mm Steel Pipe
- Side rail : 2EA
- Caster 5"



CABINET (WOOD)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|---------|
| KS-KY-056 | 400×400×700 | 1인용 상두대 |



CABINET (WOOD)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|---------------|
| KS-KY-057 | 450×420×1800 | 1인용 상두대상고 삽입형 |



CABINET (WOOD)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|------------|
| KS-KY-058 | 450×420×1800 | 개인용 Rocker |



EXAMINING BED

| Model | Size (L×D×H / mm) | Remark |
|-------|-------------------|--------|
| KS-EB | 1800×700×650 | Steel |



EXAMINING BED


| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|-------------|
| KS-EBD | 1800×700×650 | Drawer Type |



EXAMINING BED

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-EB300 | 1900×600×650 | Steel |

※ H=80T, Wood, 스펀지



가이드 레일을 내렸을때

SHOWER TROLLEY (STANDARD)

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|----------|
| KS-STC2 | 1900×600×750 | Standard |

SPECIFICATION

- ※ Stainless Steel Pipe 304 Construction.
- ※ Top plate : Long hole plate.
- ※ Side rail : swing Type 2EA.
- ※ 4inch Swivel Caster 4EA, (Brake 2EA)



INFANT BASSINET STAND

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------------|
| KS-IBS1 | 720×440×800 | High Quality |

- ※ Stainless Steel Pipe 304 Frame
- ※ P.V.C Clear crib with 50mm Mattress
- ※ 4 Inch Swivel Caster 4EA



INFANT BASSINET STAND

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|----------|
| KS-IBS2 | 720×440×800 | Standard |

- ※ Stainless Steel Pipe 304 Frame
- ※ P.V.C Clear crib with 50mm Mattress
- ※ 4 Inch Swivel Caster 4EA



STRETCHER CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|----------|
| KS-SC1 | 1900×600×700 | Standard |

- Stainless Steel 304 Construction
- 5" Swivel Brake Caster 4EA
- Accessories : Mattress, I/V Pole



STRETCHER CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-SC2 | 1900×600×700 | 동판수동형 |

- Stainless Steel 304 Construction
- 5" Swivel Brake Caster 4EA
- Accessories : Mattress, I/V Pole



SHOWER TROLLER

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-STC | 2000×750×700 | - |

FEATURE & USAGE

- 환자나 거동이 불편한 사람을 위하여 간편히 목욕 시킬 수 있다.
- 편리한 Slide One-Touch Type의 Side Rail.
- All Stainless Steel 제작으로 인한 부식의 최소화
- 방수포에 배수구 부착하여 사용상의 편리함

FUNCTIONING (Side Rail)



3단 STRETCHER CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-SC3 | 1900×700×750 | Steel |

- Body : Steel (Powder Painting Management)
- 5" Swivel Brake Caster 4EA
- Accessories : Mattress, I/V Pole



3단 STRETCHER CART

| Model | Size (L×D×H / mm) | Remark |
|--------|-------------------|--------|
| KS-SC4 | 1900×600×750 | SUS |

- Stainless Steel 304 Construction
- 5" Swivel Brake Caster 4EA
- Accessories : Mattress, I/V Pole



3단 PROTECTOR BED / CHAIR

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|------------|
| KS-YBA1003 | 561×1777×338 | BED SIZE |
| | 561×865×750 | CHAIR SIZE |



SHOWER TROLLER

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-STC1 | 2000×800×500~900 | - |

SPECIFICATION

- Stainless Steel pipe 304 Constuction, (Electronic Powder Painting Surface)
- Higher quality brake Caster 8 EA.
- Power Supply : AC 220V, 60Hz 1.5A
- Top Plate : 10mm Laminated synthetic material
- Electric code water proof Box.

FEATURE & USAGE

- 환자나 거동이 불편한 사람을 위하여 간편히 목욕 시킬 수 있으며 간병인의 편리성을 위하여 전동으로 높낮이 조절 가능.
- 전동식 리모콘을 이용하여 높이 조절이 가능하며 최저 500mm까지 내려가서 낮은 침대에도 사용이 가능하며 편리하다.
- 머리쪽에 상판의 경사각을 조절하여 샤워중이나 사용 후의 배수가 용이하다.
- 충전용 배터리와 모터는 완전 방수용으로 제작되었다.



ENDOSCOPE CART

| Model | Size (L×D×H / mm) | Remark |
|---------------|-------------------|--------|
| KS-ES300-400A | 720×600×1180 | — |
| KS-ES300-480A | 800×600×1180 | — |

※ Shelf : 400X545

SPECIFICATION

- ※ 5" Caster, Aluminum, Steel, Stainless, ABS Plastic, Rubber
- ※ Lcd Arm 적용가능
- ※ :관절형 가능, 듀얼모니터 사용 가능, 회전 및 피벗 가능
- ※ 편리한 관절식 키보드판 - 슬라이드 방식으로 대체 가능
- ※ 플라스틱 손잡이 적용 - 착용감 높임
- ※ 물림용 케이블타이 적용으로 배선 용이
- ※ 뒷문을 이용하여 장비 셋팅 및 점검 용이
- ※ 별도의 파워 스위치 사용으로 편리하고 효율적 전원관리
- ※ 중간 판넬의 상하 이동 가능
- ※ 다양한 사이즈 변경 및 제작 가능



ASSIST MONITOR CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-HDMC-01 | Φ720×400~1900 | — |

※ SPEC : 3" Caster, Wire Basket, Aluminium, Stainless, Plastic
※ 보조 모니터용 카트



ENDOSCOPE CART

| Model | Size (L×D×H / mm) | Remark |
|---------------|-------------------|---------|
| KS-ES350-400A | 560×670×1220 | 관절형 키보드 |

※ Shelf : 400X545 ※ Weight: 70kg

SPECIFICATION

- ※ 5" Caster, Aluminum, Steel, Stainless, ABS Plastic, Rubber
- ※ Lcd Arm 적용가능
- ※ :관절형 가능, 듀얼모니터 사용 가능, 회전 및 피벗 가능
- ※ 편리한 관절식 키보드판 - 슬라이드 방식으로 대체 가능
- ※ 플라스틱 손잡이 적용 - 착용감 높임
- ※ 물림용 케이블타이 적용으로 배선 용이
- ※ 뒷문을 이용하여 장비 셋팅 및 점검 용이
- ※ 별도의 파워 스위치 사용으로 편리하고 효율적 전원관리
- ※ 중간 판넬의 상하 이동 가능
- ※ 다양한 사이즈 변경 및 제작 가능



ENDOSCOPE CART

| Model | Size (L×D×H / mm) | Remark |
|---------------|-------------------|-----------|
| KS-ES350-480S | 640×670×1220 | 슬라이드형 키보드 |

※ Shelf : 480X545 ※ Weight: 70kg



- ※ 4" Caster, Aluminum, Steel, Stainless, ABS Plastic, Rubber
- ※ Load Arm 적용가능
 - : 관찰형 가능, 튜브모니터 사용 가능, 회전 및 피벗 가능
- ※ 플라스틱 손잡이 적용 - 착용감 높임
- ※ 폴링용 케이블티이 적용으로 배선 용이
- ※ 뒷문을 이용하여 장비 셋팅 및 점검 용이
- ※ 별도의 파워 스위치 사용으로 편리하고 효율적 전원관리
- ※ 중간 판매의 상하 이동 가능
- ※ 다양한 사이즈 변형 및 제작 가능



T형



7형

- Shell: 480X550 ▪ Weight: 52kg
- SPEC: 4" Caster, Keyboard Slide, 110Vx6
- ADD: Lcd Arm(s)



- * Shelf: 480X540
- * SPEC: 4" Caster, Keyboard Slide, 110V X 6
- * ADD: Lcd Arm(s)



- * Shelf: 480X545
- * SPEC: 5" Caster, Lcd Arm(s), 220V x 6, Keyboard Slide
- * ADD: IV Hnager, Oxygen tank Holder



- * Shell: 480X550
- * SPEC: 5" Caster, Lcd Arm(s), 220V x 6
110V x 6
- * ADD : Keyboard Slide, Joint Arm, Oxygen tank Holder



- Shelf:480X530
- SPEC:4 " Caster, Lcd Arm(s), 220V x 6
- ADD:IV Hnager, Keyboard Slide, Oxygen tank Holder



ARTHROSCOPE CART

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-AS-LVT-A | 520×700×1345 | - |

- * Shelf: 370X565
- * SPEC: 5" Caster, Stand Arm, Keyboard Slide
- * ADD: Joint Arm, Extension Arm, Keyboard Slide



ARTHROSCOPE CART

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-AS-LVT-B | 520×700×1345 | - |

- * Shelf: 370X565
- * SPEC: 5" Caster, Stand Arm, Keyboard Slide
- * ADD: Joint Arm, Extension Arm



일반형

램프걸이 미적용

LAPAROSCOPE CART

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-LS350-480 | 640×670×1300 | - |

- * Shelf: 480X545
- * SPEC: 5" Caster, Oxygen tank Holder, Drawer, Lcd Arm(s), Pump Hanger
- * ADD: Keyboard Slide, IV Hanger



LAPAROSCOPE CART

| Model | Size (L×D×H / mm) | Remark |
|-------------------|-------------------|--------|
| KS-LS350-480S-L0W | 640×670×1300 | - |

- * Shelf: 480X545
- * SPEC: 5" Caster, Oxygen tank Holder, Drawer, Lcd Arm(s), Keyboard Slide, Pump Hanger, IV Hanger, 220V x 6



ARTHROSCOPE CART

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-AS-AR | 550×600×1345 | - |

- * Shelf: 400X440
- * SPEC: 5" Caster, Keyboard Slide, Joint Arm, Extension Arm, Wire Basket



LAPAROSCOPE CART

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-LS300-480 | 790×600×1350 | - |

- * Shelf: 480X530
- * SPEC: 4" Caster, Oxygen tank Holder, Drawer, Lcd Arm(s), IV Hanger, 220V x 6, 110 x 6
- * ADD: Keyboard Slide, Joint Arm



LAPAROSCOPE CART

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-LS450-480 | 630×650×1255 | - |

- * Shelf: 480X530
- * SPEC: 4" Caster, Oxygen tank Holder, Drawer, Lcd Arm(s), 220V x 6
- * ADD: Keyboard Slide, IV Hanger



LAPAROSCOPE CART

| Model | Size (L×D×H / mm) | Remark |
|-----------------|-------------------|--------|
| KS-LS350-480-BB | 550×670×1300 | - |

- * Shelf: 400X545
- * SPEC: 5" Caster, Oxygen tank Holder, Drawer, Joint Arm, Extension Arm



LAPAROSCOPE CART

| Model | Size (L×D×H / mm) | Remark |
|-----------------|-------------------|--------|
| KS-LS450-480-BB | 630×670×1300 | — |

- ※ Shelf: 480X545
- ※ 4" Caster, Oxygen Tank Holder (2EA) Back



E.N.T CART

| Model | Size (L×D×H / mm) | Remark |
|--------------------|-------------------|--------|
| KS-ENT300-480-300D | 720×600×1180 | — |

- ※ Shelf: 400X530
- ※ SPEC: 5" Caster, Joint Arm, 220V×6, 110V×6
- Acrylic Hanger, Glass Tube
- ※ ADD: Drawer, Keyboard Slide



E.N.T CART

| Model | Size (L×D×H / mm) | Remark |
|----------------|-------------------|--------|
| KS-ENT300-600W | 920×600×1310 | — |

- ※ Shelf: 600X550
- ※ SPEC: 5" Caster, Joint Arm, Drawer
- ※ ADD: Lcd Arm(s), Acrylic Hanger, Glass Tube



E.N.T CART

| Model | Size (L×D×H / mm) | Remark |
|----------------|-------------------|--------|
| KS-ENT450-480W | 630×560×1045 | — |

- ※ Shelf: 480X440
- ※ SPEC: 4" Caster, 220V×6
- ※ ADD: Lcd Arm(s), Keyboard Slide, Drawer



E.N.T CART

| Model | Size (L×D×H / mm) | Remark |
|---------------|-------------------|--------|
| KS-ENT350-480 | 630×670×1090 | — |

- ※ Shelf: 480X445
- ※ SPEC: 5" Caster, Lcd Arm(s), Acrylic Hanger, 220V×6
- ※ ADD: Keyboard Slide, Drawer



E.N.T CART

| Model | Size (L×D×H / mm) | Remark |
|---------------|-------------------|--------|
| KS-ENT450-480 | 630×560×1045 | — |

- ※ Shelf: 480X445
- ※ SPEC: 4" Caster, Acrylic Hanger, 220V×6
- ※ ADD: Keyboard Slide, Drawer, Lcd Arm(s)



PC CART

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-PC350 | 680×620×835 | — |

- ※ Shelf: 520X495
- ※ SPEC: 5" Caster, Lcd Arm(s), Keyboard Slide, 220V×6
- ※ ADD: Joint Arm, Drawer



PC CART

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-PC350D | 680×620×832 | — |

- ※ Shelf: 520X495
- ※ SPEC: 5" Caster, Keyboard Slide, 220V×6
- ※ ADD: Joint Arm, Lcd Arm(s)



PC CART

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-PC450 | 670×620×800 | — |

※ Shelf : 520X495
 ※ SPEC : 4" Caster, Lcd Arm(s), Keyboard Slide, 220V×6
 ※ ADD : Joint Arm, Drawer



PC CART

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-PC350S | 680×620×1090 | — |

※ Shelf : 520X495
 ※ SPEC : 5" Caster, Keyboard Slide, 220V×6
 ※ ADD : Lcd Arm(s), Drawer Plate



PC CART

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-PC450-SDM | 670×600×1010 | — |

※ Shelf : 520X415
 ※ SPEC : 3" Caster, Lcd Arm(s), Mouse Plate, Drawer 220V x 4



PC CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-PCSU-D1 | 500×475×1045 | — |

※ 4" Caster, Lcd Arm(s), Drawer SUS 304



PC CART

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-PC350-DM | 620×570×1070 | — |

※ Shelf : 480X460
 ※ SPEC : 4" Caster, Mouse Plate, PC Bracket, 220V×6
 ※ ADD : Lcd Arm(s), Joint Arm



PC CART

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-PC450-DM | 610×520×960 | — |

※ Shelf : 470X450
 ※ SPEC : 4" Caster, Mouse Plate, PC Bracket, 220V×6
 ※ ADD : Lcd Arm(s), Joint Arm



PC CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-PCSU-DS | 500×350×1063 | — |

※ 4" Caster, Lcd Arm(s), Keyboard Slide, Drawer, SUS 304



PC CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-PC-KC01 | 510×575×885 | — |

※ Shelf : 490X460
 ※ 3" Caster, Keyboard Slide



PC CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-PC-YW01 | 600×575×880 | — |

- Shelf: 600X515
- 4" Caster, Drawer, Handle, 220V×6



PC CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-PC-KW01 | 560×580×830 | — |

- Shelf: 500X500
- 5" Caster, Handle, Keyboard Slide, 220V×6



MULTI CART

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-MT350-D2 | 590×520×815 | — |

- Shelf: 440X460
- SPEC : 4" Caster, 4 Handle, Drawer
- ADD : Keyboard Slide



MULTI CART

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-MT500-L1 | 640×520×860 | — |

- Shelf: 500X463
- SPEC : 5" Caster, 4 Handle
- ADD : Keyboard Slide, Drawer, Basket, Plate



MULTI CART

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-MT350 | 590×520×815 | — |

- Shelf: 440X460
- SPEC : 5" Caster, 4Handle, Keyboard Slide, Basket, Drawer, Plate



MULTI CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-MT350-B | 590×520×815 | — |

- Shelf: 440X460
- SPEC : 5" Caster, 4 Handle, Basket
- ADD : Keyboard Slide, Drawer



MULTI CART

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-MT500-P1 | 640×520×860 | — |

- Shelf: 500X463
- SPEC : 5" Caster, 4 Handle
- ADD : Keyboard Slide, Drawer, Basket



MULTI CART

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-MT450-D3 | 670×620×800 | — |

- Shelf: 520X495
- SPEC : 4" Caster, Drawer, 220V x 6
- ADD : Lcd Arm(s), Joint Arm



MULTI CART

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-MT450-S-B | 580×470×730 | - |

- ※ Shelf : 480×430
- ※ SPEC : 3" Caster, Keyboard Slide
- ※ ADD : Drawer, Plate



MULTI CART

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-MT350-US | 680×620×835 | - |

- ※ Shelf : 520×610
- ※ SPEC : 5" Caster, Drawer, Gel Bracket
- ※ ADD : Lcd Arm(s), Joint Arm



BOVIE CART

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-BV-BW | 610×570×1000 | - |

- ※ Shelf : 450×300, 450
- ※ 5" Caster, Drawer, Basket



BOVIE CART

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-BV-VL-UC | 500×560×880 | - |

- ※ Shelf : 470×495
- ※ 5" Caster, The Upper Shelf, Oxygen Tank Holder



MULTI CART

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-MT500-DB | 640×460×800 | - |

- ※ Shelf : 500×400
- ※ 3" Caster, Drawer, Bracket, Basket



BOVIE CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-BV450-D | 540×520×840 | - |

- ※ Shelf : 400×500
- ※ SPEC : 4" Caster, Drawer
- ※ ADD : Lcd Arm(s), Drawer Basket



BOVIE CART

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-BV-UES40 | 380×470×815 | - |

- ※ Shelf : 370×420
- ※ 4" Caster, Steel, Corner Safe Rubber



협업계용

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-SM-01 | 460×400×800 | - |

- ※ SPEC : Adjust Foot Rubber
- ※ ADD : Wire Basket, Caster



STANDARD HANGER

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|----------|
| KS-EGH-L | H=610 | Φ14 PIPE |



JOINT HANGER

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-NEGH-H | H=860 | Φ14 PIPE |

EKG CART

| Model | Cart Size (L×D×H/mm) | Inner Size (L×D×H/mm) |
|---------|----------------------|-----------------------|
| KS-EG-S | 405×325×740 | 370×285×15 |
| KS-EG-M | 530×325×740 | 490×280×15 |
| KS-EG-L | 485×405×740 | 440×365×15 |
| KS-EG-B | 540×465×740 | 490×420×15 |

※ 3" Caster, General Hanger, Steel



EKG CART (STANDERD HANGER)

| Model | Cart Size (L×D×H/mm) | Inner Size (L×D×H/mm) |
|-----------|----------------------|-----------------------|
| KS-EG-S-L | 405×325×740 | 370×285×15 |
| KS-EG-M-L | 530×325×740 | 490×280×15 |
| KS-EG-L-L | 485×405×740 | 440×365×15 |
| KS-EG-B-L | 540×465×740 | 490×420×15 |

※ 3" Caster, Standard Hanger, Steel



EKG CART (JOINT HANGER)

| Model | Cart Size (L×D×H/mm) | Inner Size (L×D×H/mm) |
|-----------|----------------------|-----------------------|
| KS-EG-S-H | 405×325×740 | 370×285×15 |
| KS-EG-M-H | 530×325×740 | 490×280×15 |
| KS-EG-L-H | 485×405×740 | 440×365×15 |
| KS-EG-B-H | 540×465×740 | 490×420×15 |

※ 3" Caster, Joint Hanger, Steel



EKG CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|----------|
| KS-EG-7402 | 485×405×820 | Hole (x) |

※ SPEC : 3" Caster, Steel
 ※ ADD : Wire Basket
 ※ 상판 : 외측 480x405x25 / 내측 440x365x15



EKG CART

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-SM-02 | 485×405×820 | - |

※ 3" Caster, Wire Basket
 ※ 상판 : 외측 480x405x25 / 내측 440x365x15



EKG CART

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|----------|
| KS-EG-7402-H | 485×405×820 | Hole (o) |

※ SPEC : 3" Caster, Steel, Paper supply space
 ※ ADD : Wire Basket
 ※ 상판 : 외측 480x405x25 / 내측 440x365x15



ECG CART

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-EG-1350 | 300×380×900 | - |

※ 4" Caster



EKG CART

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-EG-DP | 320×350×878 | - |

※ 3" Caster, Keyboard Plate, Basket



NEBULIZER

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-NE-01 | 415×450×875 | - |

※ Sun 304
※ Nebulizer



BEAUTY

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-EG450 | 310×370×733 | - |


※ 3" Caster



MEB

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-EG-MEB | 500×608×820 | - |


※ 4" Caster, Mouse Plate, Lcd Arm



ASSIST MONITOR CART (MOTOR)

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-HDMC-02 | Φ630×1430~1930 | - |


※ 2.5" Caster, Steel, Aluminium, Plastic, Motor
※ Using 220V



※ 모니터의 상하 이동을 쉽게 가능합니다.
※ 리모컨을 이용하여 높이를 조절할 수 있습니다.
※ 모니터 종류에 따라 VESA HOLE을 조절할 수 있습니다.
※ 220V의 전기를 사용합니다.



Remote control



VESA hole



Adjust hole spacing



Handle



Adjust hole spacing


모니터용 카트

환자감시용 모니터 장치 전용 카트입니다.
거의 모든 메이커를 적용할 수 있습니다.
기본적으로 형태는 같으며, 부품의 변경도 가능합니다.
체결판에 연결되는 고정판을 각 장비에 맞게 제작하여
볼트와의 체결로 부착이 됩니다.
이동과 소모품 수납이 필요한 기본 내용에 충실합니다.
수납과 물류이동 및 조립을 위해 각 부분의 분해가 가능합니다.

Air Hanger



| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-AIR-H | Φ 14 | |



인공호흡기용 카트

인공호흡기용 구성 모델입니다.
소형화를 위한 플라스틱 지름 550 미리 베이스, 안정감을 더하는 지름 120 미리 무게추, 가슴기를 부착을 위한 전용 브라켓, 각종 소모품을 수납을 위한 와이어 바구니, 두 전자기기를 이용할 수 있는 220V 멀티탭, 손쉬운 이동을 위한 손잡이가 적용됩니다. 체결판은 여러 모델들을 적용할 수 있도록 고정판과 함께 연결됩니다.



인공호흡기용 카트

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-VL-01 | Φ550×880 | - |



모니터용 카트 (고정)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-MCB-01 | Φ630×1050 | - |

높이 고정형 방식의 모델입니다. 바구니와 손잡이의 구성은 같으며, 뒷부분에는 모니터에 맞게 제작된 고정판을 체결판에 슬라이딩 하여 고정할 수 있는 마운트 판 부분이 있습니다. 모니터 고정판은 마운트 판에 삽입되어 편으로 고정됩니다. 마운트 판 부분은 약간의 틈이 가능합니다.



모니터용 카트 (고정+상하높이)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-MCB-03 | Φ630×1100~1900 | - |

높이 이동형 방식의 모델입니다. 상하 조절이 수동적으로 가능 합니다. 바구니와 손잡이의 구성은 같으며, 뒷부분에는 모니터에 맞게 제작된 고정판을 부착할 수 있는 체결판이 부착되어 있습니다. 모니터 고정판은 체결판과 함께 볼트로 체결되어 고정됩니다.



모니터용 카트 (슬라이드)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-MCB-02 | Φ630×1050 | - |

높이 고정형 방식의 모델입니다. 바구니와 손잡이의 구성은 같으며, 뒷부분에는 모니터에 맞게 제작된 고정판을 체결판에 슬라이딩 하여 고정할 수 있는 마운트 판 부분이 있습니다. 모니터 고정판은 마운트 판에 삽입되어 편으로 고정됩니다. 마운트 판 부분은 약간의 틈이 가능합니다.



모니터용 카트 (슬라이드+상하높이)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-MCB-04 | Φ630×1100~1900 | - |

높이 이동형 방식의 모델입니다. 상하 조절이 수동적으로 가능 합니다. 뒷부분에는 모니터에 맞게 제작된 고정판을 체결판에 슬라이딩 하여 고정할 수 있는 마운트 판 부분이 있습니다. 모니터 고정판은 마운트 판에 삽입되어 편으로 고정됩니다. 마운트 판 부분은 약간의 틈이 가능합니다.



체결판 (고정)
기둥의 상단에 고정됩니다. 체결판에 장비 고정판을 연결하여 사용합니다.



장비고정판 (고정)
장비 고정판은 모니터의 고정 가능한 홀의 위치에 따라 다르게 제작됩니다.



장비고정판 (슬라이드)
장비 고정판은 모니터의 특성에 따라 달라질 수 있습니다. 체결판에 삽입됩니다. 하단에 지지대가 있는 경우를 고려하여 제작되어 집니다.



01 ABS-S
작고 가볍게 사용하기에 적합합니다. 다른 베이스에 비해 중심이 조금 낮습니다.

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|---------------|
| KS-ABS-S | φ550 | A.B.S Plastic |



02 ABS-L
기본적으로 무게가 있으며 다양하게 사용됩니다. 다른 베이스에 비해 중심부분이 조금 높습니다.

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|---------------|
| KS-ABS-L | φ630 | A.B.S Plastic |



체결판 (슬라이드)
마운트 판은 좌우 삽입이 가능합니다. 핀으로 장비 부착판을 고정합니다. 마운트 판은 약간의 틸팅(각도조절)을 지원합니다.



슬라이딩 삽입
장비 고정판은 슬라이딩 방식으로 체결이 이루어집니다.



핀 고정 방식
슬라이딩 방식으로 삽입된 장비 고정판은 핀으로 고정이 이루어집니다.



03 AL-S
알루미늄 재질로 충격에 강합니다. 분체도장 방식으로 색상 변경이 가능합니다.

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|-----------|
| KS-AL-S | φ550 | Aluminium |



04 AL-L
알루미늄 재질로 충격에 강합니다. 분체도장 방식으로 색상 변경이 가능합니다.

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|-----------|
| KS-AL-L | φ700 | Aluminium |




높낮이조절 브라켓
높이를 조절할 수 있는 모델의 경우, 손잡이 부분을 돌려주어 높이를 결정한 후 잠근 후에 옆에 부착된 볼트로 2차 고정할 수 있습니다.



손잡이
이동에 편리한 손잡이는 상황에 따라 사용자에게 의해 탈부착이 가능합니다.




분리형 기둥
사용자의 이동과 물류방식에 따라 기둥을 분리하여 이동이 가능합니다. 제작방식은 사전 협의가 필요합니다.



05 무게 추
카트의 구성 후 추를 적용하여 하단 중심의 무게를 더할 수 있습니다. 안정적인 자세를 유지하기 위하여 적용합니다.

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-AC-CH120 | φ120 | 1.7kg |
| KS-AC-CH150 | φ150 | 2.6kg |
| KS-AC-CH200 | φ200 | 5kg |
| KS-AC-CH250 | φ250 | 10.5kg |



06 무게 추
무게 추는 볼트의 체결로 인해 베이스의 하단에 부착·고정됩니다.



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-IV-BHR-U | Φ550×1250~2000 | - |

※ 2.5" 바퀴, ABS베이스, 손잡이, Mash바구니, 높낮이브라켓, 보조지지대, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-IV-BHR-UC | Φ550×1250~2000 | - |

※ 2.5" 바퀴, ABS베이스, 원형손잡이+바구니, 높낮이브라켓, 보조지지대, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-IV-BHR | Φ550×1250~2000 | - |

※ 2.5" 바퀴, ABS베이스, 손잡이, Mash바구니, 높낮이브라켓, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-IV-HR | Φ550×1250~2000 | - |

※ 2.5" 바퀴, ABS베이스, 손잡이, 높낮이브라켓, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-IV-BR-U | Φ550×1250~2000 | - |

※ 2.5" 바퀴, ABS베이스, 손잡이, 높낮이브라켓, 보조지지대, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-IV-NR-U | Φ550×1250~2000 | - |

※ 2.5" 바퀴, ABS베이스, 높낮이브라켓, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-IV-BR | Φ550×1250~2000 | - |

※ 2.5" 바퀴, ABS베이스, Mash바구니, 높낮이브라켓, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-IV-NR-C | Φ550×1250~2000 | - |

※ 2.5" 바퀴, ABS베이스, 원형손잡이+바구니, 높낮이브라켓, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-IV-HR-U | Φ550×1250~2000 | - |

※ 2.5" 바퀴, ABS베이스, Mash바구니, 높낮이브라켓, 보조지지대, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-IV-NR | Φ550×1250~2000 | - |

※ 2.5" 바퀴, ABS베이스, 높낮이브라켓, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-IV-N | Φ550×1270 | - |

※ 2.5" 바퀴, ABS베이스

OPTION



선택 옵션 수액걸이대(4구) ABS조립식 타입



이동 및 조립을 위해 걸이대 부분이 접히도록 되어 있습니다.

수술실용
Operating Room
Appliances간호실용
Nursing
Appliances진료실용
Consulting Room
Appliances공급실용
Consultation Room
Appliances의료용 침대
Hospital Bed침대
Carts모니터
Monitor Carts수액 걸이대
IV Pole Stand모니터 암
Monitor Arm벽걸이형
Wall Mount폴리머
Polymer Mount액세서리
Accessory

ABS Cart

주요 제품
Order Type Carts



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|----------------|-------------------|--------|
| KS-IV-A-BH-R-U | Φ550×1250~2000 | - |

※ 2.5" 바퀴, AL베이스, 손잡이, Mash바구니, 높낮이브라켓, 보조지지대, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|----------------|-------------------|--------|
| KS-IV-A-BH-R-U | Φ550×1250~2000 | - |

※ 2.5" 바퀴, AL베이스, 원형손잡이+바구니, 높낮이브라켓, 보조지지대, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-IV-A-BH-R | Φ550×1250~2000 | - |

※ 2.5" 바퀴, AL베이스, 손잡이, Mash바구니, 높낮이브라켓, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-IV-A-H-R | Φ550×1250~2000 | - |

※ 2.5" 바퀴, AL베이스, 손잡이, 높낮이브라켓, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-IV-A-BR-U | Φ550×1250~2000 | - |

※ 2.5" 바퀴, AL베이스, 손잡이, 높낮이브라켓, 보조지지대, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-IV-A-NR-U | Φ550×1250~2000 | - |

※ 2.5" 바퀴, AL베이스, 높낮이브라켓, 보조지지대, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-IV-A-BR | Φ550×1250~2000 | - |

※ 2.5" 바퀴, AL베이스, Mash바구니, 높낮이브라켓, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-IV-A-NR-C | Φ550×1250~2000 | - |

※ 2.5" 바퀴, AL베이스, 원형손잡이+바구니, 높낮이브라켓, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|---------------|-------------------|--------|
| KS-IV-A-H-R-U | Φ550×1250~2000 | - |

※ 2.5" 바퀴, AL베이스, Mash바구니, 높낮이브라켓, 보조지지대, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-IV-A-NR | Φ550×1250~2000 | - |

※ 2.5" 바퀴, AL베이스, 높낮이브라켓, 수액걸이대(4구)



IV POLE STAND

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-IV-A-N | Φ550×1270 | - |

※ 2.5" 바퀴, AL베이스O


OPTION



선택 옵션 수액걸이대(4구) ABS조립식 타입



이동 및 조립을 위해 걸이대 부분이 접히도록 되어 있습니다.



COLUMN (S-TYPE)

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-ARM-SW | H=350 | S type |

* HOLE 19.8mm

ARM HEAD

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-ARM-HD | H=350 | Head |

- * 일반하중 : 일반 24인치 이하(포함)
- * 고 하중 : 일반 24인치 이상(포함)
- * 24인치는 공통으로 사용가능 (제조물 무게에 따라 달라질 수 있음)

모니터 거치대 (MONITOR ARM) S - TYPE

프로파일 적용 카트용 모니터 거치대(모니터암)입니다. 고산메디코 카트에 호환하여 지름은 19.8mm 홀에 삽입됩니다. VESA 규격의 홀이 있는 모니터에 대해서 거치를 지원합니다. 구매시, VESA 홀의 유무 확인이 필요합니다. 24인치 이하의 경우 360도 스윙이 가능하며, 헤드와 경우 피벗(90도)과 롤링 기능을 지원합니다.



암헤드를 기동몸체에 꽂은 후 벗어나는 것을 방지하기 위해 볼트를 조여서 빠지는 것을 방지합니다.




암헤드의 회전을 막기 위하여 노브볼트로 보조 조임 주여 회전을 방지합니다.



일반 하중
스프링의 코일 수에 따라 힘이 증가합니다. 기본으로 많이 사용되는 24인치 까지 적당합니다. 모니터가 무거우면 헤드가 내려 갑니다.



고 하중
스프링의 코일 수에 따라 힘이 증가합니다. 대화면 모니터에 사용되며 24인치 이상을 적용합니다. 모니터가 가벼우면 헤드가 올라갑니다.



모니터 거치대 I-TYPE

프로파일 적용 카트용 고산메디코 카트에 호환하여 지름 19mm에 홀에 삽입됩니다. VESA 규격의 홀이 있는 모니터에 대해서 거치를 지원합니다. 구매시, VESA 홀의 유무 확인이 필요합니다. 보다 복잡한 구조와 환경에서 거리 조정에 매우 효과적입니다. 세 부분의 센터를 중심으로 자유로운 위치 선정이 가능합니다.

EXTENSION

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-ARM-EX250 | H=250, 450 | 길이확장 |

* 단독으로 사용할 수 없고 ARM-HD와 결합 사용

COLUMN (I-TYPE)

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-ARM-I | H=300 | Column |

- * Stainless
- * ARM-EX250, 450 의 기동 역할



SIDE ARM HOLDER

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-PIH-BL | - | - |
| KS-PIH-GR | - | - |
| KS-SH100 | - | - |

- * 카트의 기동 종류의 확인이 필요합니다.
- * 300, 350, 450, 500 시리즈



고산메디코 카트에 호환하여 모니터를 카트에 추가로 거치하기 위해 필요합니다. 카트의 왼쪽 또는 오른쪽의 바깥 부분에 설치합니다. S타입과 T타입을 지원합니다. 이미 사용중인 카트에도 추가 설치 가능합니다.

JOINT ARM -A



| 1 | Model | Size (L×D×H / mm) | Remark |
|---|-------|-------------------|--------|
|---|-------|-------------------|--------|

KS-ARM-HD W=350 Head

- 일반하중 : 일반 24인치 이하(포함)
- 고 하중 : 일반 24인치 이상(포함)
- 24인치는 공통으로 사용가능
(제조품 무게에 따라 달라질 수 있음)

| 2 | Model | Size (L×D×H / mm) | Remark |
|---|-------|-------------------|--------|
|---|-------|-------------------|--------|

KS-ARM-JW W=255 -

- ARM-HD 와 결합 사용
- 관절 두 곳이 각도 90도 변화하여 높이만큼 길이의 확장이 가능합니다.

| 3 | Model | Size (L×D×H / mm) | Remark |
|---|-------|-------------------|--------|
|---|-------|-------------------|--------|

KS-ARM-I02 W=145 -

- ARM-EX 시리즈에 연결하여 기동에 결합

JOINT ARM -B



| 1 | Model | Size (L×D×H / mm) | Remark |
|---|-------|-------------------|--------|
|---|-------|-------------------|--------|

KS-ARM-HD W=350 Head

- 일반하중 : 일반 24인치 이하(포함)
- 고 하중 : 일반 24인치 이상(포함)
- 24인치는 공통으로 사용가능
(제조품 무게에 따라 달라질 수 있음)

| 2 | Model | Size (L×D×H / mm) | Remark |
|---|-------|-------------------|--------|
|---|-------|-------------------|--------|

KS-ARM-EX160 W=160 Extension

- ARM-HD 와 결합 사용

| 3 | Model | Size (L×D×H / mm) | Remark |
|---|-------|-------------------|--------|
|---|-------|-------------------|--------|

KS-ARM-JW W=255 -

- ARM-HD 와 결합 사용
- 관절 두 곳이 각도 90도 변화하여 높이만큼 길이의 확장이 가능합니다.

| 4 | Model | Size (L×D×H / mm) | Remark |
|---|-------|-------------------|--------|
|---|-------|-------------------|--------|

KS-ARM-I02 W=145 -

- ARM-EX 시리즈에 연결하여 기동에 결합

JOINT ARM -C



| 1 | Model | Size (L×D×H / mm) | Remark |
|---|-------|-------------------|--------|
|---|-------|-------------------|--------|

KS-ARM-HD W=350 Head

- 일반하중 : 일반 24인치 이하(포함)
- 고 하중 : 일반 24인치 이상(포함)
- 24인치는 공통으로 사용가능
(제조품 무게에 따라 달라질 수 있음)

| 2 | Model | Size (L×D×H / mm) | Remark |
|---|-------|-------------------|--------|
|---|-------|-------------------|--------|

KS-ARM-JW W=255 -

- ARM-HD 와 결합 사용
- 관절 두 곳이 각도 90도 변화하여 높이만큼 길이의 확장이 가능합니다.

| 3 | Model | Size (L×D×H / mm) | Remark |
|---|-------|-------------------|--------|
|---|-------|-------------------|--------|

KS-ARM-EX450 W=450 Extension

- ARM-HD 와 결합 사용
- EX250 과 결합가능

| 4 | Model | Size (L×D×H / mm) | Remark |
|---|-------|-------------------|--------|
|---|-------|-------------------|--------|

KS-ARM-I02 W=145 -

- ARM-EX 시리즈에 연결하여 기동에 결합

EXTENSION ARM -A




- 1
- ① KS-ARM-HD : 1 EA
 - ② KS-ARM-EX160 : 2 EA
 - ③ KS-ARM-I02 : 1 EA

EXTENSION ARM -B



- 2
- ① KS-ARM-HD : 1 EA
 - ② KS-ARM-EX160 : 3 EA
 - ③ KS-ARM-I02 : 1 EA

CENTER ARM -A



- 3
- ① KS-ARM-HD : 1 EA
 - ② KS-ARM-EX160 : 3 EA
 - ③ KS-ARM-C01 : 1 EA

CENTER ARM -B



- 4
- ① KS-ARM-HD : 2 EA
 - ② KS-ARM-EX160 : 6 EA
 - ③ KS-ARM-C02 : 1 EA



EXTENSION 160

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|-----------|
| KS-ARM-EX160 | W=160 | Extension |

※ ARM-HD 와 결합 사용

EXTENSION 250

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|-----------|
| KS-ARM-EX250 | W=250 | Extension |

※ ARM-HD 와 결합 사용
※ EX450과 결합가능

EXTENSION 450

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|-----------|
| KS-ARM-EX450 | W=450 | Extension |

※ ARM-HD 와 결합 사용
※ EX250과 결합가능



EXTENSION HOLDER

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-ARM-H02 | H=145 | - |

※ ARM-EX 시리즈에 연결하여 기동에 결합



CENTER HEAD (1HOLE)

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-ARM-C01 | Φ50.8 | 1 Head |

※ 모니터 ARM을 카트의 중심 부분에 거치
※ ARM-HD, EX시리즈와 연결가능



CENTER HEAD (2HOLE)

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-ARM-C02 | Φ50.8 | Head |

※ 모니터 ARM을 카트의 중심 부분에 거치
※ ARM-HD, EX시리즈와 연결가능



DUAL MONITOR ARM

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-ARM-D-1 | 문의요망 | - |



STAND MONITOR ARM

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-ARM-HDMS | 315×290×300 | - |

※ HD Medical Monitor Arm
※ 책상 거치용



[Back] 후면

TABLET HOLDER

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-T0030L | VESA100x100 | 9" ~ 13" |



[Front] 전면



[Back] 전면

TABLET HOLDER

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-T0011L | VESA100x100 | 7" ~ 11" |



NOTEBOOK PLATE (가로 슬라이드 확장)


| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-NOTE-P01 | 355x250 | - |

※ ARM-HD 와 결합 사용
※ EX250과 결합가능



SINGLE MONITOR ARM

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|---------------|
| KS-LDA-12U | VESA100x100 | LED 15" ~ 27" |





WALL MOUNT

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-WM-C | MOVE=480 | - |

- 모델별 제작 가능
- Wire Basket 추후 장착 가능

WALL MOUNT (BASKET)

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-WM-CB | MOVE=480 | — |

- 모델별 제작 가능
- Wire Basket 포함



높낮이 조절을 하여 사용자의 높이에 맞추어 사용 가능합니다.



Wire Basket은 사용자의 필요에 따라 추가 장착이 가능합니다.



Wire Basket 는 볼트의 구성만으로
추가 장착이 가능합니다.



장비의 안정된 위치를 위하여 바닥판에 장비를 고정할 수 있습니다.
회전과 약간의 틸팅을 지원하여 모델별 제작이 가능합니다.

TABLE MOUNT

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-WM-T | 바닥판 120×80 | - |

* 회전, 탈형 가능

장비 고정판은 슬라이딩 방식으로 체결이 이루어집니다.

슬라이딩 방식으로 삽입된 장비 고정판은 핀으로 고정이 이루어
집니다.



WALL MOUNT

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-WM-W | W=200 | - |

- 일반형 벽 마운트 모델
- 벽면에 LCD모니터의 고정을 하기 위하여 장착합니다.
- VESA 규격의 모니터만 가능합니다.



WALL MOUNT

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-WM-U | W=480 | - |

- 일반형 벽 마운트 모델
- 벽면에 감시용 장비의 고정을 하기 위하여 장착합니다.
- 모델별 제작이 가능합니다.



POLE MOUNT

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-PM-S | φ32 | - |
| KS-PM-L | φ51 | - |

- 원형 기둥에 추가 구성 부착 필요할 경우 사용합니다.
- 협의 제작이 필요합니다.



POLE MOUNT

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-WM-P | φ32, φ51 | VESA규격 |

- 원형 기둥에 LCD모니터를 장착하기 위해 사용합니다.



POLE MOUNT

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-HFH-38 | φ38 | - |

- For Pipe φ38
- Ventilator 에서 humidifier를 장착하기 위해 사용합니다.



PLATE MOUNT

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-PM-P | | - |

- For plate pole mount
- 플대를 사용하기 위해 테이블의 측면에 장착합니다.



WIRE BASKET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-WB-G | 300×200×150 | - |


- 모니터카트, 월마운트 등
- 용도에 따라 다양한 사용 가능



WIRE BASKET

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-WB-GS | 320×150(130)×150 | - |

- 모니터카트, 월마운트 등
- 용도에 따라 다양한 사용 가능



WIRE BASKET

| Model | Size (L×D×H / mm) | Remark |
|----------|-----------------------|--------|
| KS-WBW-W | 340(405)×210(235)×150 | - |

- 모니터카트, 월마운트 등
- 용도에 따라 다양한 사용 가능



MASH BASKET

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-MB-G | 270×150×95(140) | - |

- I.V pole stand, 모니터 카트 등
- 용도에 따라 다양한 사용 가능



WIRE BASKET

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-WB-WM | 300×200×95(150) | - |

- 다양한 용도에 사용 가능

움직임이 많고 활동이 많은 업무 가운데 필요한 임시 수납 공간을 위하여 카트 등에 Basket을 설치할 수 있습니다. 크기와 공간에 따라 적용할 수 있습니다. 기본적으로 각 Basket은 고무재질의 코팅이 되어 있습니다.



MULTI-TAP

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-220V-2 | 105×50×40 | — |

- 220V×2구, 10A
- 길이 3M



MULTI-TAP

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-220VS-2 | 130×75×40 | — |

- 220V×2구, 15A
- 길이 3M
- 개별 차단 스위치



MULTI-TAP

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-220V-4S | 235×50×40 | — |

- 220V×4구, 16A
- 길이 3M
- 통합 전원 차단 스위치



MULTI-TAP

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-220V-6S | 320×50×40 | — |

- 220V×6구, 16A
- 길이 3M (기본), 5M
- 통합 전원 차단 스위치



MULTI-TAP

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-220V-3 | 145×50×40 | — |

- 220V×3구, 10A
- 길이 3M



MULTI-TAP

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-110V-6S | 265×55×45 | — |

- 110V×6구, 16A
- 길이 3M
- 통합 전원 차단 스위치



MULTI-TAP

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-220V-10S | 490×50×40 | — |

- 220V×10구, 16A
- 길이 5M
- 통합 전원 차단 스위치



| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-IV-2 | H=1000 | - |

※ For Carts
※ $\phi 12.7, 15.8, 19.1, 25.4$

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-EH-1 | 450~550×1150 | - |

※ Endoscope Hanger
※ For 300, 350, 450, 500 Carts

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-WH-S2 | H=130 | - |

※ Endoscope Water Hanger 2 Holder



| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-CN-P | 45×45×35 | - |

※ Plate Corner Safer For Carts
※ ABS Plastic

| Model | Size (L×D×H / mm) | Remark |
|---------|-------------------|--------|
| KS-AH-B | - | - |

※ For Cart Handle
※ Steel, Aluminium, Sponge

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-GT-A1 | - | - |

※ Order Type



| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-SH-S1 | W=100 | - |

※ Endoscope Lens Safer Bracket / Holder

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-SH-S2 | W=245 | - |

※ Endoscope Lens Safer Bracket / 2 Holder

| Model | Size (L×D×H / mm) | Remark |
|-----------|--------------------|--------|
| KS-AH-TR1 | $\phi 27, \phi 38$ | - |

※ Die Casting



| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-VT-B1 | - | - |

※ Ventilator Hanger Bracket

| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-VT-B2 | 40×50×60 | - |

※ Ventilator Hanger Bracket

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-SC-125 | - | - |

※ $\phi 125$ SUS CAN
※ SUS 304



| Model | Size (L×D×H / mm) | Remark |
|--|-------------------|--------|
| KS-YDS-201 | 745x475x890 | ABS |
| ※ 선반 내측 규격(가로×세로): 550x430mm ※ 서랍 내측 규격(가로×세로×높이): 360x510x160mm ※ 최대적재량: 120kg ※ 4" 바퀴(브레이크 바퀴 x 2개) | | |



| Model | Size (L×D×H / mm) | Remark |
|--|-------------------|--------|
| KS-YDS-202 | 745x475x890 | ABS |
| ※ 선반 내측 규격(가로×세로): 550x430mm ※ 서랍 내측 규격(가로×세로×높이): 360x510x160mm ※ 최대적재량: 120kg ※ 4" 바퀴(브레이크 바퀴 x 2개) | | |



| Model | Size (L×D×H / mm) | Remark |
|--|-------------------|--------|
| KS-YDS-203 | 745x475x890 | ABS |
| ※ 선반 내측 규격(가로×세로): 550x430mm ※ 서랍 내측 규격(가로×세로×높이): 360x510x160mm ※ 최대적재량: 120kg ※ 4" 바퀴(브레이크 바퀴 x 2개) | | |



| Model | Size (L×D×H / mm) | Remark |
|---|-------------------|--------|
| KS-YDS-301 | 745x475x890 | ABS |
| ※ 선반 내측 규격(가로×세로): 550x430mm ※ 서랍 내측 규격(가로×세로×높이): 360x510x82/160mm ※ 최대적재량: 120kg ※ 4" 바퀴(브레이크 바퀴 x 2개) | | |




| Model | Size (L×D×H / mm) | Remark |
|---|-------------------|--------|
| KS-YDS-302 | 745x475x890 | ABS |
| ※ 선반 내측 규격(가로×세로): 550x430mm ※ 서랍 내측 규격(가로×세로×높이): 360x510x82/160mm ※ 최대적재량: 120kg ※ 4" 바퀴(브레이크 바퀴 x 2개) | | |



| Model | Size (L×D×H / mm) | Remark |
|---|-------------------|--------|
| KS-YDS-204 | 745x475x890 | ABS |
| ※ 선반 내측 규격(가로×세로): 550x430mm ※ 서랍 내측 규격(가로×세로×높이): 360x510x82/160mm ※ 최대적재량: 120kg ※ 4" 바퀴(브레이크 바퀴 x 2개) | | |



| Model | Size (L×D×H / mm) | Remark |
|---|-------------------|--------|
| KS-YDS-205 | 745x475x890 | ABS |
| ※ 선반 내측 규격(가로×세로): 550x430mm ※ 서랍 내측 규격(가로×세로×높이): 360x510x82/160mm ※ 최대적재량: 120kg ※ 4" 바퀴(브레이크 바퀴 x 2개) | | |



[칸막이]



| Model | Size (L×D×H / mm) | Remark |
|---|-------------------|--------|
| KS-YDS-303 | 745x475x890 | ABS |
| ※ 선반 내측 규격(가로×세로): 550x430mm ※ 서랍 내측 규격(가로×세로×높이): 360x510x82/160mm ※ 최대적재량: 120kg ※ 4" 바퀴(브레이크 바퀴 x 2개) | | |



| Model | Size (L×D×H / mm) | Remark |
|---|-------------------|--------|
| KS-YDS-304 | 745x475x890 | ABS |
| ※ 선반 내측 규격(가로×세로): 550x430mm ※ 서랍 내측 규격(가로×세로×높이): 360x510x82mm ※ 최대적재량: 120kg ※ 4" 바퀴(브레이크 바퀴 x 2개) | | |



| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-YDS-106 | 745x475x890 | ABS |

- 선반 내측 규격(가로x세로):550x430mm
- 서랍 내측 규격(가로x세로x높이):360x510x82mm
- 최대적재량:120kg
- 4" 바퀴(브레이크 바퀴 x 2개)

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-YDS-106L | 745x475x890 | ABS |

- 선반 내측 규격(가로x세로):550x430mm
- 서랍 내측 규격(가로x세로x높이):360x510x82mm
- 최대적재량:120kg
- 4" 바퀴(브레이크 바퀴 x 2개)

| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-YDS-100K | 745x475x890 | ABS |

- 선반 내측 규격(가로x세로):550x430mm
- 서랍 내측 규격(가로x세로x높이):360x510x82mm
- 최대적재량:120kg
- 4" 바퀴(브레이크 바퀴 x 2개)



| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-YDS-101 | 745x475x890 | ABS |

- 선반 내측 규격(가로x세로):550x430mm
- 서랍 내측 규격(가로x세로x높이):360x510x82mm
- 최대적재량:120kg
- 4" 바퀴(브레이크 바퀴 x 2개)

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-YDS-102 | 745x475x890 | ABS |

- 선반 내측 규격(가로x세로):550x430mm
- 서랍 내측 규격(가로x세로x높이):360x510x82mm
- 최대적재량:120kg
- 4" 바퀴(브레이크 바퀴 x 2개)

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-YDS-103 | 745x475x890 | ABS |

- 선반 내측 규격(가로x세로):550x430mm
- 서랍 내측 규격(가로x세로x높이):360x510x82mm
- 최대적재량:120kg
- 4" 바퀴(브레이크 바퀴 x 2개)



| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|----------|
| KS-YM-101A | 475x320x795 | ABS & PP |

- 선반 내측 규격(가로x세로):366x275mm
- 선반간격:555mm
- 최대적재량:100kg
- 바퀴:3" Rubber



| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|----------|
| KS-YM-101C | 475x320x1120 | ABS & PP |

- 선반 내측 규격(가로x세로):366x275mm
- 선반간격:275mm
- 최대적재량:100kg
- 바퀴:3" Rubber



| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|----------|
| KS-YM-101B | 475x320x805 | ABS & PP |

- 선반 내측 규격(가로x세로):366x275mm
- 선반간격:555mm
- 최대적재량:100kg
- 바퀴:3" Rubber



| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|----------|
| KS-YM-101D | 475x320x505 | ABS & PP |

- 선반 내측 규격(가로x세로):366x275mm
- 선반간격:275mm
- 최대적재량:100kg
- 바퀴:3" Rubber



| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|----------|
| KS-YM-101E | 475x320x435 | ABS & PP |

- 선반 내측 규격(가로x세로):366x275mm
- 선반간격:275mm
- 최대적재량:100kg



| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-101A | 560x385x820 | ABS & PP |

- 선반 내측 규격(가로x세로): 440x340mm
- 선반간격: 555mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-201A | 655x455x820 | ABS & PP |

- 선반 내측 규격(가로x세로): 525x400mm
- 선반간격: 555mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-301A | 815x490x820 | ABS & PP |

- 선반 내측 규격(가로x세로): 680x440mm
- 선반간격: 555mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-101C | 560x385x820 | ABS & PP |

- 선반 내측 규격(가로x세로): 440x340mm
- 선반간격: 475mm
- 서랍 내측 규격(가로x세로x높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-201C | 655x455x820 | ABS & PP |

- 선반 내측 규격(가로x세로): 525x400mm
- 선반간격: 475mm
- 서랍 내측 규격(가로x세로x높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-301C | 815x490x820 | ABS & PP |

- 선반 내측 규격(가로x세로): 680x440mm
- 선반간격: 450mm
- 서랍 내측 규격(가로x세로x높이): 355x260x100mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-101B | 560x385x830 | ABS & PP |

- 선반 내측 규격(가로x세로): 440x340mm
- 선반간격: 275mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-201B | 655x455x830 | ABS & PP |

- 선반 내측 규격(가로x세로): 525x400mm
- 선반간격: 275mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-301B | 815x490x830 | ABS & PP |

- 선반 내측 규격(가로x세로): 680x440mm
- 선반간격: 275mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-101D | 560x385x830 | ABS & PP |

- 선반 내측 규격(가로x세로): 440x340mm
- 선반간격: 195/275mm
- 서랍 내측 규격(가로x세로x높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-201D | 655x455x830 | ABS & PP |

- 선반 내측 규격(가로x세로): 525x400mm
- 선반간격: 195/275mm
- 서랍 내측 규격(가로x세로x높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-301D | 815x490x830 | ABS & PP |

- 선반 내측 규격(가로x세로): 680x440mm
- 선반간격: 170/275mm
- 서랍 내측 규격(가로x세로x높이): 355x260x100mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-101E | 560x385x820 | ABS & PP |

- 선반 내측 규격(가로×세로): 440x340mm
- 선반간격(가로×세로×높이): 555mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-201E | 655x455x820 | ABS & PP |

- 선반 내측 규격(가로×세로): 525x400mm
- 선반간격(가로×세로×높이): 555mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-301E | 815x490x820 | ABS & PP |

- 선반 내측 규격(가로×세로): 680x440mm
- 선반간격(가로×세로×높이): 555mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-101G | 560x385x820 | ABS & PP |

- 선반 내측 규격(가로×세로): 440x340mm
- 선반간격: 475mm
- 서랍 내측 규격(가로×세로×높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-201G | 655x455x820 | ABS & PP |

- 선반 내측 규격(가로×세로): 525x400mm
- 선반간격: 475mm
- 서랍 내측 규격(가로×세로×높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-301G | 815x490x820 | ABS & PP |

- 선반 내측 규격(가로×세로): 680x440mm
- 선반간격: 450mm
- 서랍 내측 규격(가로×세로×높이): 355x260x100mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-101F | 560x385x830 | ABS & PP |

- 선반 내측 규격(가로×세로): 440x340mm
- 선반간격(가로×세로×높이): 275mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-201F | 655x455x830 | ABS & PP |

- 선반 내측 규격(가로×세로): 525x400mm
- 선반간격(가로×세로×높이): 275mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-301F | 815x490x830 | ABS & PP |

- 선반 내측 규격(가로×세로): 680x440mm
- 선반간격(가로×세로×높이): 275mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-101H | 560x385x830 | ABS & PP |

- 선반 내측 규격(가로×세로): 440x340mm
- 선반간격: 195/275mm
- 서랍 내측 규격(가로×세로×높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-201H | 655x455x830 | ABS & PP |

- 선반 내측 규격(가로×세로): 525x400mm
- 선반간격: 195/275mm
- 서랍 내측 규격(가로×세로×높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-301H | 815x490x830 | ABS & PP |

- 선반 내측 규격(가로×세로): 680x440mm
- 선반간격: 170/275mm
- 서랍 내측 규격(가로×세로×높이): 355x260x100mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

- 브레이크 바퀴는 옵션입니다.
- 본 카탈로그의 제품사진은 실제 색상과 다소 차이가 있을 수 있으며, 예고없이 디자인이 변경 될 수 있습니다.

- 브레이크 바퀴는 옵션입니다.
- 본 카탈로그의 제품사진은 실제 색상과 다소 차이가 있을 수 있으며, 예고없이 디자인이 변경 될 수 있습니다.



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-101J | 560x385x830 | ABS & PP |

- 선반 내측 규격(가로x세로): 440x340mm
- 선반간격: 195/275mm
- 서랍 내측 규격(가로x세로x높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-201J | 655x455x830 | ABS & PP |

- 선반 내측 규격(가로x세로): 525x400mm
- 선반간격: 195/275mm
- 서랍 내측 규격(가로x세로x높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-301J | 815x490x830 | ABS & PP |

- 선반 내측 규격(가로x세로): 680x440mm
- 선반간격: 170/275mm
- 서랍 내측 규격(가로x세로x높이): 355x260x100mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-101K | 560x385x830 | ABS & PP |

- 선반 내측 규격(가로x세로): 440x340mm
- 선반간격: 195/195mm
- 서랍 내측 규격(가로x세로x높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-201K | 655x455x830 | ABS & PP |

- 선반 내측 규격(가로x세로): 525x400mm
- 선반간격: 195/195mm
- 서랍 내측 규격(가로x세로x높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-301K | 815x490x830 | ABS & PP |

- 선반 내측 규격(가로x세로): 680x440mm
- 선반간격: 170/170mm
- 서랍 내측 규격(가로x세로x높이): 355x260x100mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-401A | 560x385x930 | ABS & PP |

- 선반 내측 규격(가로x세로): 440x340mm
- 선반간격: 275/380mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-501A | 655x455x930 | ABS & PP |

- 선반 내측 규격(가로x세로): 525x400mm
- 선반간격: 275/380mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-601A | 815x490x930 | ABS & PP |

- 선반 내측 규격(가로x세로): 680x440mm
- 선반간격: 275/380mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-401B | 560x385x930 | ABS & PP |

- 선반 내측 규격(가로x세로): 440x340mm
- 선반간격: 195/380mm
- 서랍 내측 규격(가로x세로x높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-501B | 655x455x930 | ABS & PP |

- 선반 내측 규격(가로x세로): 525x400mm
- 선반간격: 195/380mm
- 서랍 내측 규격(가로x세로x높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-601B | 815x490x930 | ABS & PP |

- 선반 내측 규격(가로x세로): 680x440mm
- 선반간격: 170/380mm
- 서랍 내측 규격(가로x세로x높이): 355x260x100mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber

- 브레이크 바퀴는 옵션입니다.
- 본 카탈로그의 제품사진은 실제 색상과 다소 차이가 있을 수 있으며, 예고없이 디자인이 변경 될 수 있습니다.

- 브레이크 바퀴는 옵션입니다.
- 본 카탈로그의 제품사진은 실제 색상과 다소 차이가 있을 수 있으며, 예고없이 디자인이 변경 될 수 있습니다.



KS-Y-401C

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-401C | 560x385x930 | ABS & PP |

- 선반 내측 규격(가로×세로): 440x340mm
- 선반간격: 195/295mm
- 서랍 내측 규격(가로×세로×높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



KS-Y-501C

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-501C | 655x455x930 | ABS & PP |

- 선반 내측 규격(가로×세로): 525x400mm
- 선반간격: 195/295mm
- 서랍 내측 규격(가로×세로×높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



KS-Y-601C

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-601C | 815x490x930 | ABS & PP |

- 선반 내측 규격(가로×세로): 680x440mm
- 선반간격: 170/270mm
- 서랍 내측 규격(가로×세로×높이): 355x260x100mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



KS-Y-401E

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-401E | 560x385x930 | ABS & PP |

- 선반 내측 규격(가로×세로): 440x340mm
- 선반간격: 195/295mm
- 서랍 내측 규격(가로×세로×높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



KS-Y-501E

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-501E | 655x455x930 | ABS & PP |

- 선반 내측 규격(가로×세로): 525x400mm
- 선반간격: 195/295mm
- 서랍 내측 규격(가로×세로×높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



KS-Y-601E

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-601E | 815x490x930 | ABS & PP |

- 선반 내측 규격(가로×세로): 680x440mm
- 선반간격: 170/270mm
- 서랍 내측 규격(가로×세로×높이): 355x260x100mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



KS-Y-401D

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-401D | 560x385x930 | ABS & PP |

- 선반 내측 규격(가로×세로): 440x340mm
- 선반간격: 195/380mm
- 서랍 내측 규격(가로×세로×높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



KS-Y-501D

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-501D | 655x455x930 | ABS & PP |

- 선반 내측 규격(가로×세로): 525x400mm
- 선반간격: 195/380mm
- 서랍 내측 규격(가로×세로×높이): 355x195x75mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



KS-Y-601D

| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|----------|
| KS-Y-601D | 815x490x930 | ABS & PP |

- 선반 내측 규격(가로×세로): 680x440mm
- 선반간격: 170/380mm
- 서랍 내측 규격(가로×세로×높이): 355x260x100mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



KS-Y-701A

| Model | Size (L×D×H / mm) | Size (L×D / mm) | Remark |
|-----------|-------------------|-----------------|----------|
| KS-Y-701A | 560x385x1040 | 440x340 | ABS & PP |
| KS-Y-701B | 655x455x1040 | 525x400 | ABS & PP |
| KS-Y-701C | 815x490x1040 | 680x440 | ABS & PP |

- 선반간격: 275mm
- 최대적재량: 100kg



KS-Y-801A

| Model | Size (L×D×H / mm) | Size (L×D / mm) | Remark |
|-----------|-------------------|-----------------|----------|
| KS-Y-801A | 560x385x1120 | 440x340 | ABS & PP |
| KS-Y-801B | 655x455x1120 | 525x400 | ABS & PP |
| KS-Y-801C | 815x490x1120 | 680x440 | ABS & PP |

- 선반간격: 275mm
- 최대적재량: 100kg
- 바퀴: 4 " Rubber



| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-YS-101A | 745x475x825 | ABS |

- 선반 내측 규격(가로×세로): 550x430mm
- 선반간격: 555mm
- 최대적재량: 120kg
- 바퀴: 4" Rubber

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-YS-101B | 745x475x835 | ABS |

- 선반 내측 규격(가로×세로): 550x430mm
- 선반간격: 275mm
- 최대적재량: 120kg
- 바퀴: 4" Rubber

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-YS-101E | 745x475x825 | ABS |

- 선반 내측 규격(가로×세로): 550x430mm
- 선반간격: 275mm
- 최대적재량: 120kg
- 바퀴: 4" Rubber



| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-YS-101B-H | 745x475x835 | ABS |

- 선반 내측 규격(가로×세로): 550x430mm
- 선반간격: 275mm
- 최대적재량: 120kg
- 바퀴: 4" Rubber

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-YS-101B-W | 745x475x835 | ABS |

- 선반 내측 규격(가로×세로): 550x430mm
- 선반간격: 275mm
- 최대적재량: 120kg
- 바퀴: 4" Rubber

| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-YS-101B-C | 745x475x835 | ABS |

- 선반 내측 규격(가로×세로): 543x425mm
- 선반간격: 275mm
- 최대적재량: 120kg
- 바퀴: 4" Rubber

● 브레이크 바퀴는 옵션입니다.



| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-YS-101F | 745x475x835 | ABS |

- 선반 내측 규격(가로×세로): 550x430mm
- 선반간격: 275mm
- 최대적재량: 120kg
- 바퀴: 4" Rubber

| Model | Size (L×D×H / mm) | Size (L×D / mm) | Remark |
|----------|-------------------|-----------------|----------|
| KS-Y-100 | 560x385x535 | 440x340 | ABS & PP |
| KS-Y-200 | 655x455x535 | 525x400 | ABS & PP |
| KS-Y-300 | 815x490x535 | 680x440 | ABS & PP |

- 선반간격: 275mm
- 최대적재량: 100kg
- 바퀴: 4" Rubber

| Model | Size (L×D×H / mm) | Size (L×D / mm) | Remark |
|-----------|-------------------|-----------------|----------|
| KS-Y-100B | 560x385x435 | 440x340 | ABS & PP |
| KS-Y-200B | 655x455x435 | 525x400 | ABS & PP |
| KS-Y-300B | 815x490x435 | 680x440 | ABS & PP |

- 선반간격: 275mm
- 최대적재량: 100kg

OPTIONAL ACCESSORIES



다용도 행거(PART-SB)

IV-스탠드

행거

핸들 (HANDLE)

Y-101시리즈

Y-201시리즈

Y-301시리즈

다용도 통

거치형 쓰레기통

폐기물 거치대

형자백 수거함

내측 규격 (L×D×H) 375x120x130mm

칸막이

악게이스

● 브레이크 바퀴는 옵션입니다.



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-001 | 205×355×812 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-002 | 217×385×712 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-003 | 590×555×805 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-004 | 280×280×860 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-005 | 210×295×740 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-006 | 300×400×995 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-007 | 160×230×696 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-008 | 200×244×935 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-009 | 460×919×782 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-010 | 340×420×894 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-011 | 520×525×750 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-012 | 480×550×785 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-013 | 480×550×785 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-014 | 480×550×785 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-015 | 440×495×800 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-016 | 410×463×745 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-017 | 520×550×875 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-018 | 400×500×813 | - |

수술실용
Operating Room
Appliances간호실용
Nursing
Appliances진료실용
Consulting Room
Appliances공급실용
Supply Room
Appliances의료용 침대
Hospital Bed카트
Carts모니터 카트
Monitor Carts수액 관리대
IV Pole Stand모니터 암
Monitor Arm벽면걸
Wall Mount폴대걸
Pole Mount액세서리
AccessoryABS 카트
ABS Cart수술실용
Operating Room
Appliances



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-019 | 550×490×950 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-020 | 440×463×1250 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-021 | 520×550×800 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-028 | 480×520×1220 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-029 | 460×470×590 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-030 | 480×500×1220 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-022 | 400×400×800 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-023 | 545×430×785 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-024 | 440×500×800 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-031 | 520×500×800 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-032 | 345×395×790 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-033 | 600×520×800 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-025 | 580×550×850 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-026 | 450×500×820 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-027 | 480×623×1221 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-034 | 260×300×1133 | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-035 | — | — |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-036 | 400×400×870 | — |

수술실용
Operating Room
Appliances간호실용
Nursing
Appliances진료실용
Consulting Room
Appliances공급실용
Supply Room
Appliances의료용 침대
Hospital Bed카트
Carts모니터 카트
Monitor Carts수액 걸이대
IV Pole Stand모니터 받침
Monitor Arm벽 부착용
Wall Mount물 부착용
Pole Mount액세서리
AccessoryABS 카트
ABS Cart주요
Order Type
Carts주요
Order Type
Carts



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-037 | 480×480×870 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-038 | 480×480×870 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-039 | 480×480×870 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-046 | 460×397×950 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-047 | 433×310×800 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-048 | 480×480×870 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-040 | 420×440×1180 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-041 | 350×290×990 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-042 | 350×465×1035 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-049 | 420×480×782 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-050 | 420×480×782 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-051 | 450×520×870 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-043 | 480×390×768 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-044 | 440×380×900 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-045 | 400×540×1340 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-052 | 340×380×1064 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-053 | 430×290×800 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-054 | 920×550×1050 | - |

수술실용
Operating Room
Appliances간호
Nursing
Appliances진료실용
Consulting Room
Appliances공급실용
Supply Room
Appliances의료용 침대
Hospital Bed카트
Carts모니터 카트
Monitor Carts수액 관리대
IV Pole Stand모니터 암
Monitor Arm벽 이음돌
Wall Mount폴 이음돌
Pole Mount악세서리
AccessoryABS 카트
ABS Cart수납형 카트
Order Type Carts



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-055 | 595×525×870 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-056 | 750×800×670-970 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-057 | 320×400×855 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-064 | 410×460×1220 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-065 | 협의요망 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-066 | 협의요망 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-058 | 420×380×975 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-059 | 240×400×1050 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-060 | 360×420×800 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-067 | 450×520×873 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-068 | 300×500×745 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-069 | 920×680×765 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-061 | 520×552×768 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-062 | 440×220×750 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-063 | 800×550×1020 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-070 | 2050×770×855 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-071 | - | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-072 | 협의요망 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-073 | 900×125×900 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-074 | 1210×760×1235 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-075 | 700×440×800 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-082 | Φ550×900 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-083 | Φ434×590 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-084 | Φ550×1010 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-076 | 630×590×900 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-077 | - | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-078 | - | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-085 | Φ550×1300 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-086 | Φ550×1100 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-087 | 325×325×730 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-079 | - | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-080 | Φ550×4200~1800 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-081 | Φ550×1090 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-088 | 협의요망 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-089 | 협의요망 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-090 | 협의요망 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-091 | 홀크기 협의 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-092 | 홀위치 협의 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-093 | 협의요망 | - |



| Model | Size (L×D×H / mm) | Remark |
|----------|-------------------|--------|
| KS-NEG-H | H=860 | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-094 | - | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-095 | - | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-096 | - | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-097 | - | - |



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-098 | - | - |

수술실용
Operating Room
Appliances간호실용
Nursing
Appliances진료실용
Consulting Room
Appliances공급실용
Supply Room
Appliances의료용 침대
Hospital Bed카트
Carts모니터 카트
Monitor Carts수액 관리대
IV Pole Stand모니터 암
Monitor Arm필터문드
Wall Mount필터문드
Pole Mount악세서리
AccessoryABS 카트
ABS Cart주문형 카트
Order Type Carts



| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-099 | - | - |



| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-100 | - | - |



| Model | Size (LxDxH / mm) | Remark |
|--|-------------------|--------|
| KS-MC-CARAT | φ550×930 | - |
| * 2.5" Caster, ABS Base, ABS Handle, φ120×206 weight | | |



| Model | Size (LxDxH / mm) | Remark |
|--|-------------------|--------|
| KS-MC-DASH | φ630×910 | - |
| * 3" Caster, ABS Base, ABS Handle, Wire Basket * 모델별 상부 Plate 교체 가능 | | |



| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-101 | - | - |



| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-102 | - | 1 Head |



| Model | Size (LxDxH / mm) | Remark |
|-----------|-------------------|--------|
| KS-14-103 | - | - |



| Model | Size (LxDxH / mm) | Remark |
|--|-------------------|--------|
| KS-MC-LTV-A | φ550×1028 | - |
| * 3" Caster, ABS Base, Wire Basket, Oxygen Tank Holder | | |



| Model | Size (LxDxH / mm) | Remark |
|--|-------------------|--------|
| KS-MC-2701 | φ630×1043 | - |
| * 3" Caster, ABS Base, ABS Handle, Wire Basket | | |



| Model | Size (LxDxH / mm) | Remark |
|---|-------------------|--------|
| KS-MC-4100 | φ550×1050 | - |
| * 2.5" Caster, ABS Base, ABS Handle, Mesh Basket, φ120×20t weight | | |



| Model | Size (L×D×H / mm) | Remark |
|----------------|-------------------|--------|
| KS-MC-PM7000-H | φ630×1220~1800 | - |

※ 3" Caster, ABS Base, ABS Handle, Wire Basket

| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-MC-OXYM | φ550×1100~1600 | - |

※ 2.5" Caster, ABS Base, ABS Handle, Mesh Basket



| Model | Size (L×D×H / mm) | Remark |
|--------------|-------------------|--------|
| KS-MC-YM6000 | φ550×1092 | - |

※ 3" Caster, Die Casting Base, Wire Basket, ABS Handle



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-MC-VCM | φ550×1140~1710 | - |

※ 3" Caster, Die Casting Base, Wire Basket, ABS Handle



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-MC-ANG | φ550×1270 | - |

※ 2.5" Caster, Die Casting Base, Wire Basket
φ120×20t weight



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-MC-COD | φ700×1300 | 1 Head |

※ 3" Caster, Die Casting Base



| Model | Size (L×D×H / mm) | Remark |
|-----------|-------------------|--------|
| KS-MC-NIC | φ550×1097 | - |

※ 3" Caster, ABS Base, 220Vx2



| Model | Size (L×D×H / mm) | Remark |
|-------------|-------------------|--------|
| KS-MC-VS100 | φ550×1000~150 | - |

※ 2" Caster, ABS Base, Wire Basket, ABS Handle,
φ120×20t weight



| Model | Size (L×D×H / mm) | Remark |
|------------|-------------------|--------|
| KS-MC-T100 | φ550×930 | - |

※ 2.5" Caster, ABS Base, ABS Handle, Wire Basket,
φ120×20t weight



수술실용
Operating Room
Appliances

간호실용
Nursing
Appliances

진료실용
Consulting Room
Appliances

공급실용
Supply Room
Appliances

의료용 침대
Hospital Bed

카트
Carts

모니터 카트
Monitor Carts

수역 선반
IV Pole Stand

모니터 암
Monitor Arm

물 연결부
Water Mount

물 연결부
Pole Mount

악세서리
Accessory

ABS 카트
ABS Cart

수문용 카트
Order Type Carts




64



65



66



67




77



78



79



80



68



69



70



71



81



82



83



72



73



74



75



76



수술실용
Operating Room
Appliances
간호실용
Nursing
Appliances
진료실용
Consulting Room
Appliances
공급실용
Supply Room
Appliances
의료용 침대
Hospital Bed
가동
Carts
모니터 카트
Monitor Carts
수액ポール
IV Pole Stand
모니터 암
Monitor Arm
물리온도
Water Mount
물리온도
Pole Mount
액세서리
Accessory
ABS 카트
ABS Cart
주문형 카트
Order Type Carts



<http://www.kosanmedico.co.kr> | <http://www.kosanmedico.net>

